



ATLANTIC
MEDICAL
SOLUTIONS




Status™


COVID-19/Flu A&B

A Rapid Immunoassay for the Simultaneous
Direct Detection and Differential Diagnosis of SARS-CoV-2,
Influenza Type A and Type B Antigen from nasopharyngeal swab specimens

Infections with these viruses may present similar symptoms.
Can you tell them apart? **We Can.**



Learn More




1 test + 1 specimen = 3 Results

The additive power of **Status™ COVID-19/Flu A&B**




- **Sensitivity and Specificity**
 - COVID-19 - PPA 94%, NPA 100%
 - Flu A - PPA 89.2%, NPA 99.4%
 - Flu B - PPA 86.4%, NPA 99%
- **FDA Emergency Use Authorization (EUA)**
- **Visually read in 15 minutes**
- **Flocked nasopharyngeal swab for superior specimen collection and patient comfort**

ORDERING INFORMATION		
PRODUCT	# OF TESTS	PRODUCT CODE
Status™ Covid-19/Flu A&B	25	LSN33225

Proudly made in the USA



Easy 3 Step Procedure

<p>1</p>  <p>squeeze</p>	<p>Tear the tab off the Extraction Reagent Capsule and dispense entire contents into the Extraction Well.</p>
<p>2</p> 	<p>Insert the specimen swab in the Swab Stand.</p> <ul style="list-style-type: none"> • Rotate swab 3 times to mix the specimen. • Let stand 1 minute. • Rotate swab 3 times again and discard the swab.
<p>3</p>  <p>Tap</p>	<p>Raise the device upright and let stand 1–2 seconds.</p> <p>Gently tap the device to ensure the liquid flows into the hole.</p> <p>Lay the device back down.</p>
<p>Set a timer for 15 minutes.</p>	
<p>4</p>	<p>Read test results at 15 minutes.</p> <p>NOTE: False positive or false negative results can occur if the test is not read between 15 and 20 minutes.</p>

- Status™ COVID-19/FLU A&B is an antigen test for use with nasopharyngeal swab specimens.
- Status™ COVID-19/FLU A&B test has been granted emergency use authorization (EUA) by

the FDA for use at the point-of-care setting. During the duration of the emergency this test can be performed in a patient care setting that is operating under a CLIA Waiver, Certificate of Compliance or Certificate of Accreditation.

 **Proudly made in the USA**