






ATLANTIC  
MEDICAL  
SOLUTIONS



# PRODUCT CATALOG

got gloves?™



# We Are Your Global GLOVE EXPERT

With MORE THAN 100 YEARS of glove manufacturing experience, we have long been a pioneer within the disposable glove industry. We are able to leverage the latest manufacturing technology in **research** and **innovation** to ensure our gloves achieve optimal **safety**, **reliability**, and **quality**.



## WORLD CLASS MANUFACTURING

**Setting standards of excellence.**

Our glove facilities are some of the most sophisticated manufacturing plants in the world. We utilize state-of-the-art technology to achieve the highest levels of quality in our medical and industrial grade gloves.



## RESEARCH & DEVELOPMENT

**Science never stands still.**

Neither do we. We are constantly conducting advanced research and development in our ongoing commitment to develop new products that will complement our existing lines.



## COMPLIANCE & INTEGRITY

**We are honest and fair.**

We handle our resources carefully and treat our customers, employees, suppliers, and environment with respect. We are focused on safety and strive to continuously improve working conditions to retain our highly qualified work force.



Rev01122024

## PEOPLE, PROTECTING PEOPLE

At Sempermed, we are proud helping people. Proud of our role as protectors. And proud to now be part of HARPS Global.



[www.gotgloves.com](http://www.gotgloves.com) • Ph 1.800.366.9545



# CHOOSING THE RIGHT GLOVE

GLOVE MATERIAL	ATTRIBUTES	ISSUES	PROPERTIES	APPLICATIONS
<b>Chloroprene</b> <i>For maximum protection, change gloves every 30 minutes</i> 	<ul style="list-style-type: none"> <li>Excellent barrier protection</li> <li>Greater puncture resistance than latex</li> <li>Comfort: similar to latex</li> <li>Finish: textured</li> </ul>	<ul style="list-style-type: none"> <li>Accelerators in glove manufacturing process may cause skin irritation (Less common compared to latex allergies)</li> </ul>	<ul style="list-style-type: none"> <li>Molecular structure allows wearer to quickly identify a rip, hole or tear</li> <li>Displays good flame resistance and performs well in a wide range of temperatures</li> </ul>	IDEAL FOR PROLONGED USE & APPROPRIATE FOR MOST APPLICATIONS medical, dental, extreme temperatures, chemical exposure, lab, dermatology, jan/san, safety, and beauty
<b>Nitrile</b> <i>For maximum protection, change gloves every 20 - 30 minutes</i> 	<ul style="list-style-type: none"> <li>Excellent barrier protection</li> <li>Greater puncture resistance than latex</li> <li>Comfort: similar to latex</li> <li>Available in a variety of colors</li> <li>Finish: textured</li> </ul>	<ul style="list-style-type: none"> <li>Accelerators in glove manufacturing process may cause skin irritation (Less common compared to latex allergies)</li> </ul>	<ul style="list-style-type: none"> <li>Molecular structure allows wearer to quickly identify a rip, hole or tear</li> <li>Less susceptible to material breakdown and leakage than a vinyl or latex glove</li> </ul>	APPROPRIATE FOR MOST APPLICATIONS medical, dental, vet, chemical exposure, lab, jan/san, safety, food services, beauty, and tattoo
<b>Latex</b> <i>For maximum protection, change gloves every 20 - 30 minutes</i> 	<ul style="list-style-type: none"> <li>Excellent barrier performance</li> <li>Appropriate for most uses</li> <li>Highly elastic, comfortable fit</li> <li>Finish: micro-textured and/or fully textured</li> </ul>	<ul style="list-style-type: none"> <li>Potential for allergies to latex proteins. Reaction can range from sneezing to anaphylactic shock (life-threatening)</li> <li>Allergies affect 10 to 20% of work-exposed individuals</li> </ul>	<ul style="list-style-type: none"> <li>Strength and elasticity</li> <li>Barrier protection against a variety of chemicals</li> <li>Less susceptible to material breakdown and leakage</li> </ul>	APPROPRIATE FOR MOST APPLICATIONS general patient examinations, dental, wound care, routine oral procedures, lab, food services, beauty, and jan/san
<b>Vinyl</b> <i>For maximum protection, change gloves every 15 minutes or less</i> 	<ul style="list-style-type: none"> <li>Most economical choice</li> <li>Comfort: not critical</li> <li>Available in a variety of colors</li> <li>Finish: smooth</li> </ul>	<ul style="list-style-type: none"> <li>Cuff has loose fit and is susceptible to material breakdown and leakage</li> <li>Natural body oils can breakdown the glove material</li> </ul>	<ul style="list-style-type: none"> <li>Minimal elasticity/flexibility</li> <li>Typically clear in color, however colored vinyl is widely used in food manufacturing</li> </ul>	INTENDED FOR SHORT-TERM NON-INVASIVE TASKS dispensing medication, non-invasive physical examinations, jan/san, food services, and food manufacturing
<b>Polyethylene</b> <i>For maximum protection, change gloves frequently</i> 	<ul style="list-style-type: none"> <li>Most economical choice for non-medical markets</li> <li>Finish: smooth</li> </ul>	<ul style="list-style-type: none"> <li>Baggy, loose fit</li> <li>Not intended for medical use or stringent chemicals</li> <li>Very thin material</li> <li>Intended for very short-term tasks only</li> </ul>	<ul style="list-style-type: none"> <li>Minimal elasticity/flexibility</li> <li>Typically clear in color</li> <li>Low protection</li> <li>Susceptible to ripping/tearing and leakage</li> </ul>	INTENDED FOR LIGHT DUTY APPLICATIONS jan/san, food services, cleaning * NOT FOR MEDICAL USE



SUSA REORDER #	PRODUCT NAME	GRADE	CHEMO TESTED	ESD TESTED	SKIN FRIENDLY	FENTANYL	CBD OIL	GASTRIC ACID	SPECIAL COATING	GLOVE COLOR
NITRILE										
BKNF	SemperForce® Black	Exam	—	✓	—	✓	✓	—	—	Black
BLNF	SemperForce® Blue	Exam	✓	✓	—	✓	—	—	—	Blue
GSBNF	GripStrong® Black	Industrial	—	—	—	—	—	—	—	Black
GSBLN	GripStrong® Blue	Industrial	—	✓	—	—	—	—	—	Blue
GSWNF	GripStrong® White	Industrial	—	✓	—	—	—	—	—	White
INIHD	SemperGuard® HD	Industrial	—	✓	—	—	—	—	—	Blue
INIPFT	SemperGuard®	Industrial	—	✓	—	—	—	—	—	Blue
NIPFT	SemperCare®	Exam	✓	✓	—	✓	—	—	—	Blue
ORNF	SemperForce® Orange	Exam	—	—	—	✓	—	—	—	Orange
SMNP	StarMed® PLUS	Exam	✓	—	—	✓	—	—	—	Purple
SMNR	Best Touch® StarMed® Rose	Exam	—	—	—	—	—	—	Aloe	Pink
SMNS	StarMed® Select	Exam	✓	—	—	✓	—	—	—	Purple
SMTN	StarMed® Ultra	Exam	✓	—	—	✓	—	—	—	Purple
SSEC	SemperShield® Extended Cuff	Exam	✓	—	—	✓	—	—	—	Blue
SSSC	SemperShield® Standard Cuff	Exam	—	✓	—	✓	—	—	—	Blue
SUNB	SemperSure® Black	Exam	✓	✓	✓	✓	—	✓	Accelerator-Free	Black
SUNF	SemperSure®	Exam	✓	—	✓	✓	—	✓	Accelerator-Free	Blue
SUNG	SemperSure®	Exam	✓	—	✓	—	—	—	Accelerator-Free	Green
TTNF200	TenderTouch® 200	Exam	✓	—	—	✓	✓	—	—	Blue
TTNF100	TenderTouch® 100	Exam	✓	—	—	✓	✓	—	—	Blue
LATEX										
GSLF	GripStrong®	Industrial	—	—	—	—	—	—	—	Orange
INDPFT	SemperGuard®	Industrial	—	✓	—	—	—	—	—	Orange
INDPS	SemperGuard®	Industrial	—	—	—	—	—	—	Powdered	Orange
SCLT	SemperCare®	Exam	—	—	—	—	—	—	—	Orange
SM	StarMed®	Exam	—	—	—	—	—	—	—	Orange
VINYL										
EVNP	SemperCare® Synthetic	Exam	—	—	—	—	—	—	—	White
SCVNP	SemperCare® Stretch	Exam	—	—	—	—	—	—	—	Orange
GSVF	GripStrong®	Industrial	—	—	—	—	—	—	—	White
VBLK	SemperGuard® Black	Industrial	—	—	—	—	—	—	—	Black
VP	SemperGuard®	Industrial	—	—	—	—	—	—	Powdered	White
VPF	SemperGuard®	Industrial	—	—	—	—	—	—	—	White
POLYETHYLENE										
HPEF	SemperGuard® FoodSafe Stretch Poly	Industrial	—	—	—	—	—	—	—	White



# Glove Chemical Resistance Chart\*

\*This chart is provided as a guideline only. It does not suggest that any Sempermed® glove has been tested for use with the chemical listed. Sempermed® accepts no responsibility for end-user glove selection against a specific risk. Conditions of use and appropriate glove choice are the responsibility of the user.

E – Excellent    G – Good    F – Fair    P – Poor    N – Not Rated



CHEMICAL	NITRILE (NBR)	NATURAL RUBBER LATEX (NRL)	VINYL (PVC)	CHLOROPRENE (NEOPRENE)
Acetaldehyde	P	G	G	G
Acetic Acid	F	E	E	E
Acetone	P	G	P	P
Alcohol (Isopyl)	E	E	E	E
Alcohol (Ethyl)	E	E	G	E
Alcohol (Isopropyl)	E	E	G	E
Ammonium Hydroxide (conc.)	E	G	E	F
Aniline	E	F	F	G
Ammonia	P	G	P	P
Antifreeze (methanol-based)	E	E	P	E
Benzaldehyde	F	F	G	F
Benzene	P	P	F	F
Benzyl Chloride	N	F	P	P
Bleach	P	G	P	P
Brake Cleaner (containing ethanol or hexane)	E	P	P	P
Brake Cleaner, non-chlorinated (containing acetone, n-heptane or xylene)	E	P	P	P
Bromine	N	G	G	G
Calcium Hypochlorite	G	P	G	G
Carbon Disulfide	G	P	P	P
Carbon Tetrachloride	F	P	F	F
CBD Oil	E	N	N	N
Chlorine	E	G	G	G
Chloroacetone	N	F	P	E
Chloroform	P	P	P	F
Chromic Acid	G	P	E	F
Cyclohexane	E	F	P	E
Cytotoxic Chemicals	E	G	P	E
Dibenzyl Ether	N	F	P	G
Dibutyl Phthalate	G	F	P	G
Diesel Fuel	E	P	P	E
Diethanolamine	G	F	E	E
Diethyl Ether	E	F	P	G
Dimethyl Sulfoxide	E	N	N	N
Ethyl Acetate	P	F	P	G
Ethylene Dichloride	P	P	P	F
Ethylene Glycol	E	G	E	G
Fentanyl	E	G	P	E
Fluorine	N	G	G	G
Formaldehyde	E	G	E	E
Formic Acid	E	G	E	E
Gasoline (oil)	E	P	P	E

CHEMICAL	NITRILE (NBR)	NATURAL RUBBER LATEX (NRL)	VINYL (PVC)	CHLOROPRENE (NEOPRENE)
Gastric Acid	G	P	P	G
Grease	G	P	P	G
Glucose	G	P	P	G
Glycerol	E	G	E	G
Hair Dye	G	P	P	N
Hand Sanitizer	E	E	G	E
Hexane	E	P	P	E
Hydrochloric Acid (conc.)	G	P	F	F
Hydrofluoric Acid (50%)	F	P	F	E
Hydrogen Peroxide	G	G	E	G
Insecticides / Pesticides	P	G	F	P
Iodine	G	G	G	E
Kerosene	E	P	P	E
Methylamine	E	G	E	G
Methyl Cellosolve	P	F	P	E
Methyl Chloride	P	P	P	E
Methylene Chloride	P	F	F	F
Methyl Ethyl Ketone	P	F	P	G
Mineral Spirits	P	G	F	P
Monoethanolamine	N	F	E	E
Morpholine	P	F	E	E
Motor Oil	G	P	P	N
Naphthalene	E	G	G	G
Nitric Acid (conc.)	P	P	G	P
Ortho-phthalaldehyde	G	P	P	G
Paint (water based)	G	P	P	N
Paint (oil-based)	G	P	P	N
Perchloric Acid	F	P	E	G
Peroxide	G	P	P	N
Phenol	F	G	E	G
Phosphoric Acid	E	G	E	E
Potassium Hydroxide (wet)	E	G	E	G
Propane	E	P	P	N
Sodium Hydroxide	E	G	E	E
Sodium Hypochlorite	E	G	G	E
Sulfuric Acid (conc.)	F	G	G	E
Toluene	P	P	F	F
Trichloroethylene	P	P	F	F
Tricresyl Phosphate	E	P	P	F
Triethanolamine	E	F	E	E
Turpentine	E	P	P	N
Xylene	F	P	N	P

REV 04/202022



# NITRILE Disposable Gloves



## SemperGuard<sup>®</sup> HD Nitrile

### FEATURES

- Industrial grade heavy duty nitrile
- ESD tested and approved
- Fully textured
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 6-7 mil glove



### Order # Size

INIHD102	S
INIHD103	M
INIHD104	L
INIHD105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

72 (12 x 6)

All Sempermed gloves are single use only.



## SemperGuard<sup>®</sup> Nitrile

### FEATURES

- Industrial grade blue nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 4-5 mil glove



### Order # Size

INIPFT102	S
INIPFT103	M
INIPFT104	L
INIPFT105	XL
INIPFT106	XXL

### GLOVE COUNT

100 Gloves/Box (90 XXL)  
10 Boxes/Case

### PALLET CONFIGURATION

78 (13 x 6)

All Sempermed gloves are single use only.

REV01162024



Industrial Grade



ESD Tested & Approved



Food Handling Compliant

sempermed

semperguard

100+ YEARS  
- EST. 1920 -





# NITRILE Disposable Gloves



## GRIPSTRONG® BLUE NITRILE

### FEATURES

- Industrial grade blue nitrile
- Textured fingertips
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

GSBLN102 S

GSBLN103 M

GSBLN104 L

GSBLN105 XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.



## GRIPSTRONG® WHITE NITRILE

### FEATURES

- Industrial grade white nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

GSWNF102 S

GSWNF103 M

GSWNF104 L

GSWNF105 XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.



## GRIPSTRONG® BLACK NITRILE

### FEATURES

- Industrial grade black nitrile
- Textured fingertips
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

GSBNF102 S

GSBNF103 M

GSBNF104 L

GSBNF105 XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

78 (13 x 6)

All Sempermed gloves are single use only.

REV01162024



Industrial Grade



ESD Tested & Approved



Food Handling Compliant





# NITRILE Disposable Gloves



## SEMPERSHIELD®

### FEATURES

- Exam grade nitrile - Fully textured
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 6-7 mil glove

*For more information about protection against fentanyl, please contact your sales rep.*



### Order # Size

SSSC102	S
SSSC103	M
SSSC104	L
SSSC105	XL
SSSC106	XXL

### GLOVE COUNT

100 Gloves/Box (90 XXL)  
10 Boxes/Case

### PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



## SEMPERSHIELD® EXTENDED CUFF

### FEATURES

- Exam grade nitrile - Extended cuff
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Tested for use against Fentanyl
- **Thickness:** Competes with a 5-6 mil glove

*For more information about protection against fentanyl, please contact your sales rep.*



### Order # Size

SSEC102	S
SSEC103	M
SSEC104	L
SSEC105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

78 (13 x 6)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



ESD Tested & Approved



Extended Cuff



Chemo Tested



Food Handling Compliant







# NITRILE Disposable Gloves



## SEMPERFORCE **BLACK**

### FEATURES

- Exam grade black nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

BKNF102	S
BKNF103	M
BKNF104	L
BKNF105	XL
BKNF106	XXL

### GLOVE COUNT

100 Gloves/Box (90 XXL)  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



## SEMPERFORCE **BLUE**

### FEATURES

- Exam grade blue nitrile
- Textured fingertips
- Meets USP 800 requirements for PPE compliance for chemotherapy
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

BLNF102	S
BLNF103	M
BLNF104	L
BLNF105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



## SEMPERFORCE **ORANGE**

### FEATURES

- Exam grade orange nitrile
- Fully textured
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

ORNF102	S
ORNF103	M
ORNF104	L
ORNF105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



ESD Tested & Approved



Chemo Tested



CBD Oil



Food Handling Compliant





# NITRILE Disposable Gloves



## StarMed® ULTRA Nitrile Exam

### FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- **Thickness:** Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

SMTN251 XS

SMTN252 S

SMTN253 M

SMTN254 L

SMTN255 XL

### GLOVE COUNT

250 Gloves/Box (225 XL)  
10 Boxes/Case

### PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



## StarMed® PLUS Nitrile Exam

### FEATURES

- Exam grade nitrile
- Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- **Thickness:** Competes with a 3 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

SMNP301 XS

SMNP302 S

SMNP303 M

SMNP304 L

SMNP305 XL

### GLOVE COUNT

300 Gloves/Box (250 XL)  
10 Boxes/Case

### PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



## StarMed® Select Nitrile Exam

### FEATURES

- Exam grade nitrile
- Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- **Thickness:** Competes with a 3 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

SMNS101 XS

SMNS102 S

SMNS103 M

SMNS104 L

SMNS105 XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

64 (16 x 4)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



ESD Tested & Approved



Chemo Tested



Food Handling Compliant







# NITRILE Disposable Gloves



Best Touch®  
**StarMed® ROSE**  
Nitrile Exam with Aloe

## FEATURES

- Exam grade nitrile
- Textured fingertips
- Coated with aloe
- Moisturizes dry, chapped hands
- **Thickness:** Competes with a 3-4 mil glove



## Order # Size

SMNR201 XS

SMNR202 S

SMNR203 M

SMNR204 L

SMNR205 XL

## GLOVE COUNT

200 Gloves/Box (180 XL)  
10 Boxes/Case

## PALLET CONFIGURATION

72 (12 x 6)

All Sempermed gloves are single use only.



## Tender Touch® 200

Nitrile Exam

## FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



## Order # Size

TTNF201 XS

TTNF202 S

TTNF203 M

TTNF204 L

TTNF205 XL

## GLOVE COUNT

200 Gloves/Box (180 XL)  
10 Boxes/Case

## PALLET CONFIGURATION

72 (12 x 6)

All Sempermed gloves are single use only.



## Tender Touch® 100

Nitrile Exam

## FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



## Order # Size

TTNF101 XS

TTNF102 S

TTNF103 M

TTNF104 L

TTNF105 XL

## GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

## PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



Chemo Tested



CBD Oil



Food Handling Compliant



Aloe Coated

sempermed

semperguard

100+ YEARS  
EST. 1920



# NITRILE Disposable Gloves



## SemperSure® Black

### FEATURES

- Accelerator-Free exam grade nitrile
- Tested for protection against Fentanyl & gastric acid
- Meets USP 800 requirements for PPE compliance for chemotherapy
- ESD tested & approved
- **Thickness:** Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

SUNB101	XS
SUNB102	S
SUNB103	M
SUNB104	L
SUNB105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

64 (16 x 4)

All Sempermed gloves are single use only.



## SemperSure®

### FEATURES

- Accelerator-Free exam grade nitrile
- Textured fingertips
- Low dermatitis potential
- Tested for protection against Fentanyl & gastric acid
- Meets USP 800 requirements for PPE compliance for chemotherapy
- **Thickness:** Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



### Order # Size

SUNF201	XS
SUNF202	S
SUNF203	M
SUNF204	L
SUNF205	XL

### GLOVE COUNT

200 Gloves/Box (180 XL)  
10 Boxes/Case

### PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



## SemperSure® Green

### FEATURES

- Accelerator-Free exam grade nitrile
- Textured fingertips
- Innovative green technology
- No chlorine used in manufacturing process
- Meets USP 800 requirements for PPE compliance for chemotherapy
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

SUNG201	XS
SUNG202	S
SUNG203	M
SUNG204	L
SUNG205	XL

### GLOVE COUNT

200 Gloves/Box (180 XL)  
10 Boxes/Case

### PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.

REV01162024



# 5 REASONS to Choose SemperSure® Green

## 2 MANUFACTURED WITHOUT CHLORINE

- Removing chlorine in the mfg. process is more eco-friendly, as this process creates harmful by-products that can be released into the air and water.
- Chlorination can cause gloves to become more brittle over time, reducing their elasticity. It can also leave a residual odor.\*



## 3 EARTH-FRIENDLY GREEN MANUFACTURING

- Sustainably produced with our innovative, green technology conserving natural resources:
  - 12% less water consumption
  - 17% less CO<sub>2</sub> emissions
- Bright, mint-green color to emphasize the eco-friendly benefits.



## 4 ENHANCED COMFORT

- Advanced nitrile formulation facilitates optimal 'touch and feel' instilling confidence when handling medical instruments.
- Promotes easy donning thanks to an inner polymer coating.



## 1 SKIN-FRIENDLY ACCELERATOR-FREE

- Manufactured without harmful chemical accelerators - reduces risk of skin irritation.
- Low dermatitis claim recognized by FDA and tested on humans.
- Preventative care to help reduce the risk of developing sensitivity to nitrile gloves (contact dermatitis) due to prolonged glove use.



## 5 CHEMO CLAIM

- Meets USP 800 requirements for PPE compliance for chemotherapy.
- To carry a chemo claim, the glove must adhere to specific requirements set forth by the FDA and pass testing\* against permeation of the actual cytotoxic chemicals.



LET'S SAVE THE WORLD TOGETHER



Our inner boxes are made of 80% recycled material.  
Our outer cartons are made of 90% recycled material.

\* Testing results available by request.



Glove the Planet® is our commitment to reducing our environmental impact. By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit [www.glovetheplanet.com](http://www.glovetheplanet.com).





# LATEX Disposable Gloves



## SemperCARE®

### FEATURES

- Exam grade latex
- Micro-rough textured surface
- **Thickness:** Competes with a 5-6 mil glove



\*Not applicable where prohibited.

### Order # Size

SCLT101	XS
SCLT102	S
SCLT103	M
SCLT104	L
SCLT105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



Aloe Coated



Food Handling Compliant

sempermed

semperguard

100+ YEARS  
- EST. 1920 -





# LATEX Disposable Gloves



## SemperGuard® Latex Powder-Free

### FEATURES

- Industrial grade latex
- Micro-rough textured surface
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 5-6 mil glove



\*Not applicable where prohibited.

Order #	Size
---------	------

INDPFT101	XS
-----------	----

INDPFT102	S
-----------	---

INDPFT103	M
-----------	---

INDPFT104	L
-----------	---

INDPFT105	XL
-----------	----

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



## SemperGuard® Latex Powdered

### FEATURES

- Industrial grade latex
- Powdered
- Micro-rough textured surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 4-5 mil glove



\*Not applicable where prohibited.

Order #	Size
---------	------

INDPS102	S
----------	---

INDPS103	M
----------	---

INDPS104	L
----------	---

INDPS105	XL
----------	----

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

112 (16 x 7)

All Sempermed gloves are single use only.

REV01162024



Industrial  
Grade



ESD Tested  
& Approved



Powdered



Food Handling  
Compliant

semper  
med

semper  
guard

100+  
YEARS  
- EST. 1920 -



# VINYL Disposable Gloves



## SemperCARE® VINYL STRETCH

### FEATURES

- Exam grade vinyl
- Smooth surface
- Stretch properties and cream color
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 4-5 mil glove



### Order # Size

SCVNP101	XS
SCVNP102	S
SCVNP103	M
SCVNP104	L
SCVNP105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



## SemperCARE® VINYL Synthetic

### FEATURES

- Exam grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

EVNP101	XS
EVNP102	S
EVNP103	M
EVNP104	L
EVNP105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



## GRIPSTRONG® VINYL

### FEATURES

- Industrial grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



### Order # Size

GSVF102	S
GSVF103	M
GSVF104	L
GSVF105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

112 (16 x 7)

All Sempermed gloves are single use only.

REV01162024



Exam Grade



Industrial Grade



Food Handling Compliant







# VINYL Disposable Gloves



## SemperGuard® Vinyl Powder Free

### FEATURES

- Industrial grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



Order #	Size
VPF102	S
VPF103	M
VPF104	L
VPF105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

96 (16 x 6)

All Sempermed gloves are single use only.



## SemperGuard® Vinyl Powdered

### FEATURES

- Industrial grade vinyl
- Smooth surface
- Powdered
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3-4 mil glove



Order #	Size
VP102	S
VP103	M
VP104	L
VP105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

112 (16 x 7)

All Sempermed gloves are single use only.



## SemperGuard® Black Vinyl

### FEATURES

- Industrial grade black vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 3 mil glove



Order #	Size
VLK102	S
VLK103	M
VLK104	L
VLK105	XL

### GLOVE COUNT

100 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

112 (16 x 7)

All Sempermed gloves are single use only.

REV01162024



Industrial Grade



Powdered



Food Handling Compliant

sempermed

semperguard

100+ YEARS  
- EST. 1920 -



# POLYETHYLENE Disposable Gloves



## SemperGuard® FoodSafe

### FEATURES

- Industrial grade polyethylene
- Stretch properties
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- **Thickness:** Competes with a 2 mil glove



Order #	Size
HPEF202	S
HPEF203	M
HPEF204	L
HPEF205	XL
HPEF206	XXL

### GLOVE COUNT

200 Gloves/Box  
10 Boxes/Case

### PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.

REV01162024



Industrial  
Grade



Food Handling  
Compliant

semper  
med

semper  
guard

100+  
YEARS  
- EST. 1920 -



# SURGICAL Disposable Gloves

## sempermed® supreme



### GLOVE COUNT

50 Pairs/Box ; 6 Boxes/Case

### FEATURES

- Sterile Surgical Glove
- Latex
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65

• **Thickness:** Competes with a 4-5 mil glove  
*Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.*

### PALLET CONFIG.

32 (8 x 4)



Order #	Size
SPFP550	5.5
SPFP600	6
SPFP650	6.5
SPFP700	7
SPFP750	7.5
SPFP800	8
SPFP850	8.5
SPFP900	9

All Sempermed gloves are single use only.



## sempermed® syntegra cr



### GLOVE COUNT

40 Pairs/Box ; 6 Boxes/Case

### PALLET CONFIG.

36 (16 x 6)

### FEATURES

- Sterile Surgical Glove
- Chloroprene
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65
- **Thickness:** Competes with a 5-6 mil glove

Order #	Size
SCR600	6
SCR650	6.5
SCR700	7
SCR750	7.5
SCR800	8
SCR850	8.5

All Sempermed gloves are single use only.

## sempermed® syntegra IR



### GLOVE COUNT

50 Pairs/Box ; 6 Boxes/Case

### PALLET CONFIG.

32 (8 x 4)



Order #	Size
SIR550	5.5
SIR600	6
SIR650	6.5
SIR700	7
SIR750	7.5
SIR800	8
SIR850	8.5
SIR900	9

All Sempermed gloves are single use only.



### FEATURES

- Sterile Surgical Glove
- Polyisoprene
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65
- **Thickness:** Competes with a 4-5 mil glove

REV01162024



# 100+ YEARS

- EST. 1920 -

PEOPLE, PROTECTING PEOPLE

REV01162024



**got gloves?**<sup>™</sup>

[sempermedusa.com](http://sempermedusa.com) | [gotgloves.com](http://gotgloves.com)  
1-800-366-9545



# WE DELIVER

## People, Protecting People



### DEDICATED GLOVE EXPERT

Gloves are our ONLY business. We dipped our first glove 100 years ago. Sempermed has expertise in every phase of the manufacturing process. The result is a product you can trust from people who know their business.

### RELIABLE SALES SUPPORT

You can count on your dedicated Account Manager to be there. We are committed to help resolve any issues and find solutions for you and your customers.



### DOMESTIC LOGISTICS TEAM




Sempermed's logistics team works with our factories and transportation vendors daily, so you don't have to. SUSA operates in the EST time zone.

### CREATIVE MARKETING SUPPORT



Need help with your marketing or branding? Not a problem! Our talented marketing gurus can help with private label artwork and custom marketing pieces tailored to fit your needs. Artwork fees apply.



### NO LANGUAGE BARRIER

We communicate with our factories on a daily basis speaking factory "language" for efficient order processing.



### IN HOUSE PRODUCT SOURCING

Looking for a specific glove? No problem. Let us find it for you.



### DOOR TO DOOR DELIVERY

Saving working capital while a shipment is in transit, terms begin when it lands at your door.

### HASSLE FREE FDA & CUSTOMS CLEARANCE



Relax. Let our team of experts get your gloves across the border.

### WEEKLY UPDATES

We offer weekly updates on every order you place to keep you informed.



## DEDICATED TO QUALITY PRODUCTS, SERVICE & SUPPORT

At Sempermed, we manufacture disposable gloves from start to finish. With more than 100 years of knowledge and experience, we are one of the most efficient and innovative glove manufacturers in the world. We control the process from R&D, production, and logistics to the end user. We are experts at packaging and delivering our full line of products including latex, nitrile, vinyl, chloroprene, and polyethylene gloves. Few companies can promise that kind of dedication to your hand protection. That's why we say Sempermed gloves will meet your expectations, or they will never reach your hands. We are your global glove expert, one that can be trusted unconditionally.



People, Protecting People.  
[www.gotgloves.com](http://www.gotgloves.com)  
Ph 1.800.366.9545





# HANDHYGIENE

## Infection Prevention

The average  
**HEALTH CARE  
WORKER HAND**  
contains up to  
**4.6 MILLION**  
CFUs  
(colony forming units)


**80%**

of illness  
causing  
germs are  
**SPREAD BY  
YOUR HANDS**



### GLOVES ARE **NOT** ENOUGH

The World Health Organization (WHO) recommends properly washing hands **BEFORE & AFTER GLOVE USE** by rubbing hands with alcohol-based hand sanitizer or by washing with soap and water



### WHEN SHOULD YOU CHANGE YOUR GLOVES?

- 1 As soon as gloves are **DAMAGED**
- 2 **AFTER CONTACT** with a patient or customer and his/her surroundings
- 3 **AFTER CONTACT** with blood, saliva, or other bodily fluids
- 4 When there is evidence of **CONTAMINATION**



### PREVENT GERM TRANSMISSION & GLOVE CROSS CONTAMINATION

PROPERLY  
don & doff  
gloves



LIMIT the #  
of surfaces  
you touch  
with gloved  
hands

DISPOSE  
gloves  
safely



### CHANGE GLOVES OFTEN

### HOW OFTEN SHOULD YOU CHANGE YOUR GLOVES?

For maximum protection, it's recommended to change your gloves every:



**CHLOROPRENE**  
20 - 30 MINS



**NITRILE**  
30 MINS



**LATEX**  
15 - 20 MINS



**VINYL**  
15 MINS  
OR LESS



**POLYETHYLENE**  
5-10 MINS  
OR LESS

On average,  
you come in  
contact with  
**300  
SURFACES**  
every  
**30 MINUTES**,  
exposing  
you to  
**840,000  
GERMS**

### MINIMIZE THE EXPOSURE. MINIMIZE THE RISK.

Let's face it, the world's a pretty germey place. So we've made it our sole mission to provide our customers protection from all those wicked whatevers. Since we focus only on manufacturing and supplying gloves, we set the bar on quality standards that surpass regulatory guidelines for both consistency and performance. Because there's a person attached to every single one of the hands we protect - we feel that it's people who we're really protecting. We take enormous pride in being the first and last line of defense for people in harm's way.

Source: World Health Organization, 2009.

**S semper med**

**S semper guard**



**100+ YEARS**  
EST. 1920


People. Protecting People.  
[www.gotgloves.com](http://www.gotgloves.com)  
Ph 1.800.366.9545





# FACING FENTANYL

Stay Informed. Stay Protected.



FENTANYL WAS INVENTED IN **1960** for treating severe pain due to illnesses such as cancer or as a surgical anesthetic



Fentanyl is often **ACCIDENTALLY ABUSED** & commonly mixed with heroin or cocaine in **POWDER, PILL, OR LIQUID FORM** and is often manufactured in China.


Fentanyl is a **HIGHLY ADDICTIVE** synthetic opioid that is **80 TO 100 TIMES** more potent than morphine



FENTANYL IS SO POTENT that a dose the size of a **GRAIN OF SAND** can be deadly



YOU CAN'T SEE, SMELL OR TASTE FENTANYL



Fentanyl can easily be **ABSORBED THROUGH THE SKIN**, putting medical staff, first responders, police, security, coroners, and morticians **AT RISK**




For the first time in US history, our nation surpassed **100,000 OVERDOSE FATALITIES** in a 12 month period.

**OVER 64,000** of those deaths **WERE DUE TO FENTANYL.**



Fentanyl has quickly & quietly **BECOME THE LEADING CAUSE OF DEATH** among adults **AGES 18-45.**



The CDC confirms that Fentanyl is now the **DEADLIEST DRUG IN AMERICA**

**80%** of Fentanyl in America comes into our country **THROUGH THE US-MEXICO BORDER**



For more information about protection against fentanyl, please contact your sales rep.



## got FENTANYL gloves?

As a global glove expert, it is our job to educate about the risks you may face on the job, such as Fentanyl exposure. To help minimize the risk, Sempermed offers a handful of nitrile exam gloves that provide protection against harmful fentanyl contact. Choose from SemperShield® Standard Cuff, SemperShield® Extended Cuff, SemperForce® Orange, SemperForce® Blue, SemperForce® Black, StarMed® PLUS, StarMed® Ultra, StarMed® Select, Tender Touch® 100, Tender Touch® 200, SemperSure® and SemperSure® Black.

Testing for protection against Fentanyl is conducted by an independent lab and verified by SUSA Regulatory. Testing results are available by request.

Sources: <https://www.murphy.house.gov>, <https://www.dea.gov/fentanylawareness>



People. Protecting People.  
[www.gotgloves.com](http://www.gotgloves.com)  
Ph 1.800.366.9545





# CHEMOGLOVES

Stay Informed. Stay Protected.

## PEOPLE, PROTECTING PEOPLE.



According to the World Health Organization, there were **14.1 MILLION** **NEW CANCER CASES & 32.6 MILLION PEOPLE LIVING WITH CANCER** (within 5 years of diagnosis) in 2012 worldwide



At present, more than **100** **ANTICANCER DRUGS** are used in daily clinical practice



Our Chemotherapy tested gloves meet **USP 800** **REQUIREMENTS FOR PPE COMPLIANCE FOR CHEMOTHERAPY**



In order for a glove to carry a chemo claim, the glove must adhere to specific requirements set forth by the **FDA** & PASS TESTING AGAINST **PERMEATION** of the actual cytotoxic chemicals



The **EXCRETIONS** OF PATIENTS such as urine, feces, vomit, and sweat, may be **CONTAMINATED** with chemo drugs, making urine bags, clothing and hospital beds hazardous objects



## got chemo gloves?

We have made it our mission to provide our customers with protection from all those wicked whatevers. Since there is a person attached to every single one of the hands we protect, we feel it is the people themselves that we are really protecting. If you handle cytotoxic chemicals or work in oncology, try one of our nitrile exam glove brands that have been tested for protection against chemotherapy drugs\*: SemperShield® Extended Cuff, SemperSure® Black, SemperSure® Green, SemperSure® SemperForce Blue®, TenderTouch® 200, TenderTouch® 100, StarMed® PLUS, StarMed® Ultra, and StarMed® Select.

\*Meets USP 800 requirements for PPE compliance for chemotherapy

**S** semper  
med

**S** semper  
guard

**HARPS**  
GLOBAL

**100** YEARS  
EST. 1910

People, Protecting People.  
[www.goigloves.com](http://www.goigloves.com)  
Ph 1.800.366.9545





# ESD GLOVES

We've Got You Covered



## WHAT IS ESD

### ELECTROSTATIC DISCHARGE:

The sudden flow of electricity between two electrically charged objects caused by contact, an electrical short, or dielectric breakdown




The human body is an **EXCELLENT CONDUCTOR OF ELECTRICITY** and a person can become statically charged with upwards of **SEVERAL THOUSAND VOLTS**



The static charge that builds up in humans must be released, usually resulting in **AN ELECTRIC SHOCK**

By definition, a static dissipative material has a **SURFACE RESISTIVITY** of  $1 \times 10^5$  to  $1 \times 10^{12}$  ohms per square



In many industries, this discharge of static electricity can **DAMAGE A VARIETY OF PRODUCTS**, especially sensitive electronics

## STATIC ELECTRICITY CAN POSE A FIRE-RISK



in industries that handle flammable liquids, such as automotive, painting, chemical, lab, pyrotechnic and petroleum/gas

## ESD GLOVES

provide **HAND PROTECTION** that **HELPS PREVENT PRODUCT DAMAGE** from electrostatic discharge and are **GENERALLY SAFE** to work with flammable liquids/chemicals or static sensitive devices



# got ESD gloves?

Are you handling sensitive electronics or flammable liquids? We've got you covered! Sempermed offers ten brands that are ESD tested and approved. ESD gloves provide hand protection that helps prevent product damage from electrostatic discharge, and they are safe to work with static sensitive devices where bare hands are not desired\*. We can send you free samples of our ESD tested and approved gloves: SemperSure® Black Nitrile, SemperForce® Black Nitrile, SemperForce® Blue Nitrile, SemperShield® standard cuff, SemperGuard® Nitrile, SemperGuard® HD Nitrile, GripStrong® Blue Nitrile, GripStrong® White Nitrile, and SemperGuard® PF Latex.

\* Level of ESD protection should be confirmed for your needs.

semper  
med

semper  
guard



100  
YEARS

People. Protecting People.  
www.gotgloves.com  
Ph 1.800.366.9545






# SENSITIVE SKIN?

We've Got You Covered




**SENSITIVITY TO LATEX PROTEINS**  
remains a primary cause of concern among glove users

**LATEX PROTEIN ALLERGIES AFFECT 4.3%**  
of the general population




**ALOE**  
has a soothing effect on the skin and has been shown to **HELP REDUCE ITCHINESS & INFLAMMATION**

**TYPE I (IMMEDIATE)**  
Hypersensitivity to natural rubber latex proteins can be **VERY SERIOUS OR OCCASIONALLY EVEN FATAL**



**TYPE IV (DELAYED)**  
Hypersensitivity which involves a skin rash and is often **CAUSED BY CHEMICALS** used in the manufacturing process




**SOME GLOVE USERS** can develop rashes linked to irritant contact **DERMATITIS**. This may be caused by chemicals such as **ACCELERATORS** used in the manufacture of the gloves




Many hospitals and healthcare facilities have **MIGRATED AWAY FROM LATEX GLOVES**, yet some users experience skin sensitivities to non-latex gloves

**TRY ACCELERATOR-FREE GLOVES TO MINIMIZE THE RISK**



## ARE YOU SENSITIVE TO YOUR GLOVES?

Do you have delicate skin or suffer from dry, chapped hands? Do you often wonder if you are sensitive to your gloves? We've got you covered! Sempermed offers a handful of glove brands that provide relief for irritated hands. Best Touch® StarMed® ROSE nitrile exam has an aloe coating, which soothe and nourish skin irritation. Our SemperSure® glove family are accelerator-free, low-dermatitis, nitrile exam gloves designed for those with skin sensitivities to standard nitrile gloves. They offer all the benefits of a nitrile glove, minus the chemical accelerators (added in the manufacturing process) that may cause skin reactions.

COATED WITH ALOE



ACCELERATOR-FREE



ACCELERATOR-FREE



ACCELERATOR-FREE



**S semper med**

**S semper guard**



**100+ YEARS**  
EST. 1920

People. Protecting People.  
[www.gotgloves.com](http://www.gotgloves.com)  
Ph 1.800.366.9545



Rev02232024

Skin-Friendly Gloves



# GIVING BACK

Sempermed Proudly Supports



## Clothes To Kids



### \$65 = 1 WARDROBE

Since 2003, Sempermed has donated more than

**4,054**  
WARDROBES

#### A WARDROBE INCLUDES:

- 5 pairs of new underwear
- 5 pairs of new socks
- 5 tops (including uniforms)
- 4 bottoms (including uniforms)
- 1 dress (optional)
- 1 pair of shoes
- 1 jacket (in season)

### GRASS ROOTS

3 store locations in Florida  
& affiliate locations in  
CO, NC, CT, & RI



THE HIGHEST RATING  
given to a charity

**89%**

OF EXPENDITURES  
go directly to the program

### IT ALL STARTED WITH ONE PURPLE COAT

Read the full story by visiting  
[www.clothestokids.org](http://www.clothestokids.org)

Clothes To Kids envisions a community in which **EVERY SCHOOL-AGE CHILD HAS QUALITY CLOTHING** so that they may attend school with the **CONFIDENCE & SELF-ESTEEM** needed to achieve academic success



## 2024 ANNUAL CHARITY EVENTS



APRIL 4  
5:30PM-9PM

### 22nd Birthday - Let's Shop

Armature Works  
An evening of shopping, food, open bar,  
live music, and fundraising



JULY & AUGUST  
(SUSA COLLECTS IN MAY)

Drop Off Your Drawers  
New Underwear Drive



OCTOBER

Shoetober - New Shoe Drive &  
Swing Fore Shoes Event

## CLOTHE A KID. CHANGE A LIFE.

Clothes To Kids was founded in 2002 in Pinellas County, Florida, home of Sempermed USA's headquarters. Their mission is to provide new & used quality clothing to low-income or in crisis school-age children, free of charge. \$65 will provide one child with a wardrobe; that's a full week's worth of school clothes. **Want to help? Visit [www.clothestokids.org](http://www.clothestokids.org). All donations are tax-deductable.**



People. Protecting People.  
[www.gotgloves.com](http://www.gotgloves.com)  
Ph 1.800.366.9545





# SUSTAINABILITY

## Contributing to a Prosperous world

### COMPLIANCE & INTEGRITY

are the foundation of our sustainability strategy  
SEMPERMED'S  
5 CORE PILLARS  
OF SUSTAINABILITY:

1. RESOURCE MGMT AND ENVIRONMENT
2. HEALTH AND SAFETY
3. EMPLOYEES AND SOCIETY
4. SUPPLIERS
5. INNOVATION



### SAFETY FIRST

We strive to PREVENT ACCIDENTS and continuously IMPROVE WORKING CONDITIONS to retain our highly qualified, enthusiastic workforce. To achieve this, we remain focused on **TRANSPARENCY, INTEGRITY AND CONSISTENCY**



Business Social Compliance Initiative

The Semperit Group has committed itself to high ethical and social standards, as well as fair working conditions for all employees world-wide. We are a proud BSCI member and proponent of employee rights.

The objective of the EcoVadis Corporate Social Responsibility (CSR) rating is to measure the quality of a company's CSR management system through its policies, actions, and results.



2022: Semperit Group has achieved the **PLATINUM STATUS ECOVADIS CSR RATING** - Top 1% of all companies rated



2021: Semperit Group has achieved the **GOLD STATUS ECOVADIS CSR RATING**.



2020 & 2019: Our Wimpassing site was awarded the **GOLD STATUS ECOVADIS CSR RATING** for performance in the area of corporate social responsibility.



**HEALTH, SAFETY, ENVIRONMENT, ENERGY and QUALITY (HSEQ)** are integral parts of our business.

### THIS SEAL IS OUR PROMISE FOR:

- Sustainably produced skin-friendly gloves
- Strict quality controls for optimal protection
- Reliability by ensuring European standards
- Use of renewable energies in the production process
- Innovative mfg. processes and products



### YOU CAN'T MANAGE WHAT YOU DON'T MEASURE


Semperit participates in the international CARBON DISCLOSURE PROJECT (CDP) emissions reporting program, joining over 9,600 other companies who report on CLIMATE CHANGE, SUPPLY CHAIN, WATER SECURITY & FORESTS.



We focus on encouraging all employees to achieve continuous improvement and are committed to our long-term goal of **ZERO INCIDENTS**. We strive to **AVOID WASTE** and higher consumption in order to ensure **RESOURCE-EFFICIENT, SUSTAINABLE MANUFACTURING OF FIRST-CLASS PRODUCTS**.

## COMMITTED TO SUSTAINABLE GROWTH

We are aware of our responsibility towards customers, employees, suppliers, shareholders, society, and the environment. We handle our resources carefully, treat our fellow human beings with respect, and consider the consequences of our actions at all times. That is why we have anchored sustainability within our strategic corporate goals and formulated correspondingly strict standards. Our group-wide sustainability policy is based on our code of conduct. This is supplemented by globally applicable guidelines and standards in the five focus areas - resource management, health and safety, work and society, suppliers, and innovation. We are taking clear steps towards a sustainable company development. We can make a difference together.



got gloves?

People. Protecting People.  
www.gotgloves.com  
Ph 1.800.366.9545

Follow Us!



SCAN ME



# That's how we **GLOVE** *the* **PLANET**®

GLOVE THE PLANET® IS OUR COMMITMENT TO REDUCING OUR ENVIRONMENTAL IMPACT.

By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations.



## Recycled Packaging

### We Save Trees!

- We use over **5 million pounds of recycled materials** in our packaging every year.
- Our recycled packaging **saves approximately 5,000 trees** annually.



## Reduce Waste

### Taking Steps in the Right Direction

- Our goal is to **reduce our material waste by 4% each year.**
- We produce flip flops from **leftover nitrile glove scraps.**



## Green Manufacturing

### We Save Resources

- Using no chlorine and **12% less** water and gas in our innovative green mfg. process, reduces our **CO<sub>2</sub> impact by more than 17%.**



## 30 by 2030

### We Save Energy

- We have a mission to **reduce** our overall energy, emissions, waste, and water **by 30%** by the year 2030. That's a **3% reduction each year.**



**got gloves?**™

LET'S **SAVE**  
**THE WORLD**  
TOGETHER



**GLOVE** *the* **PLANET**®

SCAN ME



SCAN TO LEARN MORE ABOUT  
OUR COMMITMENT TO SUSTAINABILITY

Rev07052023

**sempermed**®  
YOUR GLOBAL GLOVE EXPERT™

**100** YEARS  
EST. 1920

Our Gloves are Disposable. The People they Protect are Not.®

[www.sempermedusa.com](http://www.sempermedusa.com) • Ph 1.800.366.9545



**GLOVE the PLANET®**

[www.GlovethePlanet.com](http://www.GlovethePlanet.com)

## OUR COMMITMENT TO REDUCING OUR ENVIRONMENTAL IMPACT

By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations.

### RECYCLABLE

Alternative to conventional waste disposal by converting uncontaminated waste materials into new materials.

Non-contaminated recycled gloves can be up-cycled into new eco-friendly plastic products.



### COMPOSTABLE

Compostable plastics are those that biodegrade in industrial composting operations at the rate of other compostable materials.

Compostable gloves eventually break down into dirt.



### BIODEGRADABLE

All organic matter will eventually biodegrade. The term biodegradable must be qualified by the environment and timeframe.

Biodegradable gloves eventually break down into small pieces.

Biodegradable does not mean a material is compostable or recyclable.



Only Non-contaminated gloves can be recycled or sent to the land fill. Gloves in contact with blood, body fluids or chemicals must be disposed of as waste.



**LATEX GLOVES**

Can take  
1-5 Years\*  
to break down



**NITRILE GLOVES**

Can take  
Decades\*  
to break down



**VINYL GLOVES**

Can take  
over 100 years\*  
to break down



**BIODEGRADABLE GLOVES**

Can take  
up to 5 years\*  
to break down

\*Depends on landfill conditions such as humidity & temperature



### RAW MATERIALS ARE LIMITED

Since 1970 the global raw material consumption has more than tripled.

### ENERGY IS VALUABLE

By 2050 the global energy demand will increase by 70%-100%.

### CLIMATE CHANGE

Is thought to be the #1 threat to global public health today.

Sources: Definitions: [www.merriam-webster.com](http://www.merriam-webster.com), [www.businessdictionary.com](http://www.businessdictionary.com) Department of Ecology (WA), Focus on "Biobased," "Biodegradable," & "Compostable" Plastics, Aug. 2014

**sempermed®**  
YOUR GLOBAL GLOVE EXPERT™

**100+ YEARS**  
EST. 1920

Our Gloves are Disposable. The People they Protect are Not.®  
Ph 1.800.366.9545 | [www.sempermedusa.com](http://www.sempermedusa.com)  
[www.glovetheplanet.com](http://www.glovetheplanet.com)

LET'S **SAVE**  
**THE WORLD**  
TOGETHER

SCAN ME




**gotgloves®**

REV02062023



# got gloves?



## QUALITY GLOVES BY **sempermed®** THAT LEAVE A **LASTING IMPRESSION**

ENHANCED COMFORT & FIT • SUPERIOR PROTECTION • LESS HAND FATIGUE



### **SEMPERFORCE® BLACK** Nitrile Exam

REORDER# **BKNF102-106**

Sizes **S - XXL**

- Strong & durable
- Excellent gripping capabilities

100 gloves/bx (90 XXL) | 10 bxs/cs



### **StarMed® ROSE** Nitrile Exam with Aloe

REORDER# **SMNR201-205**

SIZES **S - XL**

- Coated w/ Aloe to hydrate dry hands
- Helps keep cuticles soft

200 gloves/bx (180 XL) | 10 bxs/cs



### **SemperSure® Black** Nitrile Exam

REORDER# **SUNB102-105**

SIZES **S - XL**

- Accelerator-free/low-derm
- Ideal for sensitive skin

100 gloves/bx | 10 bxs/cs



GLOVE the PLANET®  
www.GloveThePlanet.com

Glove the Planet® is our commitment to reducing our environmental impact. By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit [www.glovetheplanet.com](http://www.glovetheplanet.com).

PRODUCT  
CATALOG



SCAN ME

REV0304/2024



[AtlanticMedicalSolutions.com](http://AtlanticMedicalSolutions.com)

704-561-0004