



# T-Cup<sup>®</sup> Multi-Drug Urine Test Cup



- 2-26 configurations available
- Results in 5 minutes
- Easy to use
- 24 month shelf life
- 3 adulterations offered: CR, SG, PH
- More adulterations optional for OEM products

## CLIA Waived Drugs & Cutoff Levels

AMP AMP AMP AMP AMP AMP AMP AMP	Amphetamine (AMP)	1000ng/mL 500ng/mL	TCA TCA TCA TCA TCA TCA TCA TCA	Tricyclic Antidepressants (TCA)	1000ng/mL
BAR BAR BAR BAR BAR BAR BAR BAR	Barbiturates (BAR)	300ng/mL	MET MET MET MET MET MET MET MET	Methamphetamine (MET/mAMP)	1000ng/mL 500ng/mL
BZO BZO BZO BZO BZO BZO BZO BZO	Benzodiazepines (BZO)	300ng/mL	MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA	Methylenedioxymethamphetamine (MDMA)	500ng/mL
BUP BUP BUP BUP BUP BUP BUP BUP	Buprenorphine (BUP)	10ng/mL	MOP MOP MOP MOP MOP MOP MOP MOP	Morphine (MOP/OPI300)	300ng/mL
COC COC COC COC COC COC COC COC	Cocaine (COC)	300ng/mL 150ng/mL	OPI OPI OPI OPI OPI OPI OPI OPI	Opiate (OPI)	2000ng/mL
THC THC THC THC THC THC THC THC	Cannabinoids (THC)	50ng/mL	OXY OXY OXY OXY OXY OXY OXY OXY	Oxycodone (OXY)	100ng/mL
EDDP EDDP EDDP EDDP EDDP EDDP EDDP EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300ng/mL	PCP PCP PCP PCP PCP PCP PCP PCP	Phencyclidine (PCP)	25ng/mL
MTD MTD MTD MTD MTD MTD MTD MTD	Methadone (MTD)	300ng/mL	PPX PPX PPX PPX PPX PPX PPX PPX	Propoxyphene (PPX)	300ng/mL

## Forensic Use Drugs & Cutoff Levels

6-MAM 6-MAM 6-MAM 6-MAM 6-MAM 6-MAM 6-MAM 6-MAM	Monoacetylmorphine (6-MAM)	10ng/mL	K2 K2 K2 K2 K2 K2 K2 K2	Synthetic Cannabinoid (K2)	50ng/mL
AMP AMP AMP AMP AMP AMP AMP AMP	Amphetamine (AMP)	300ng/mL 500ng/mL	KET KET KET KET KET KET KET KET	Ketamine (KET)	300ng/mL 1000ng/mL
BZO BZO BZO BZO BZO BZO BZO BZO	Benzodiazepines (BZO)	100ng/mL 200ng/mL	KRA KRA KRA KRA KRA KRA KRA KRA	Kratom (KRA)	300ng/mL
COC COC COC COC COC COC COC COC	Cocaine (COC)	100ng/mL 150ng/mL	MET MET MET MET MET MET MET MET	Methamphetamine (MET/mAMP)	300ng/mL 500ng/mL
COT COT COT COT COT COT COT COT	Cotinine (COT)	200ng/mL	MOP MOP MOP MOP MOP MOP MOP MOP	Morphine (MOP/OPI)	100ng/mL
EDDP EDDP EDDP EDDP EDDP EDDP EDDP EDDP	EDDP	100ng/mL	MTD MTD MTD MTD MTD MTD MTD MTD	Methadone (MTD)	200ng/mL
EtG EtG EtG EtG EtG EtG EtG EtG	Ethyl Glucuronide (EtG)	500ng/mL	MQL MQL MQL MQL MQL MQL MQL MQL	Methaqualone (MQL)	300ng/mL
FTY FTY FTY FTY FTY FTY FTY FTY	Fentanyl (FTY)	20ng/mL	THC THC THC THC THC THC THC THC	Cannabinoids (THC)	25ng/mL 40ng/mL
HMO HMO HMO HMO HMO HMO HMO HMO	Hydromorphone (HMO)	300ng/mL	TRA TRA TRA TRA TRA TRA TRA TRA	Tramadol (TRA)	100ng/mL 200ng/mL 1000ng/mL

Generic stock product adulterates include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.

## Ordering Information

※ Box Size / Case Size: 25 / 4

Item Code	Description	Complexity
TDOA-154	5-Panel T-Cup; AMP, COC, OPI, PCP, THC	CLIA WAIVED
TDOA-254	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-254A3	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC with 3 AD	CLIA WAIVED
TDOA-164	6-Panel T-Cup; AMP, COC, MET, OPI, PCP, THC	CLIA WAIVED
TDOA-264	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-274	7-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
TDOA-3104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA WAIVED
TDOA-3104A3	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with 3 AD	CLIA WAIVED
TDOA-4104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-7104	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA WAIVED
TDOA-8104	10-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA WAIVED
TDOA-3124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-4124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, PCP, PPX, THC	CLIA WAIVED
TDOA-6124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-6124A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-6125A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125V	12-Panel T-Cup w/Vacu Lid; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-7125	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-7125A3	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-2135	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-2135A3	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-6135F	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA	Forensic Use Only
TDOA-7135A3EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, ETG, FTY, TRA with 3 AD	Forensic Use Only
TDOA-8135E	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA WAIVED
TDOA-1144A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-1145A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-8145EF	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only
TDOA-9145EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2	Forensic Use Only