

NitraCare Nitrile Exam Gloves

Powder-Free Non-Sterile Nitrile Examination Gloves

NitraCare Nitrile Exam Gloves are our best-selling exam gloves in our full Exam glove portfolio. They provide outstanding barrier protection from bacteria, viruses and chemicals during procedures where risk of fluid exposure is moderate to high.

Characteristics and Overview:

- Finger-tip textured enable grip even when wet.
- Standard color coded sizing to identify different sizes.
- Offer excellent sensitivity for touch perception and dependable protection.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs*



Item No	Size	Color	Packaging
MGLMG505XS	X-Small	Cobalt Blue	200 gloves/Box 2000 gloves/Case
MGLMG505S	Small	Cobalt Blue	200 gloves/Box 2000 gloves/Case
MGLMG505M	Medium	Cobalt Blue	200 gloves/Box 2000 gloves/Case
MGLMG505L	Large	Cobalt Blue	200 gloves/Box 2000 gloves/Case
MGLMG505XL	X-Large	Cobalt Blue	200 gloves/Box 2000 gloves/Case
MGLMG505XXL	2X-Large	Cobalt Blue	180 gloves/Box 1800 gloves/Case



Glove Specifications:

Brand	NitraCare
Material	Nitrile
Color	Cobalt Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	XS, SM, MD, LG, XL and 2XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Textured Fingers
Length (inches)	9.5 in
Palm Thickness (Avg)	3.0 mil
Finger Thickness (Avg)	3.6 mil
Tensile Strength (Mpa)	35 / 33 Mpa
Ultimate Elongation (%)	600 / 560
Pin Hole (AQL)	1.5



Quality Standards:

Manufactured in accordance with Quality System ISO 9001.
 Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978.
 This product is made from Synthetic Nitrile Polymer and contains no natural latex rubber.
 * When used in accordance to the guidelines set forth in the USP 800 chapter.
 Tested chemotherapy drugs available upon request.



NitraCare Nitrile Exam Gloves - MG505

Powder-Free Non-Sterile Nitrile Examination Gloves

A comfortably soft feel with an outstanding tactile sensitivity.

Chemo & Fentanyl Tested



Meets ASTM D6978 - 05

Standard Practice for Assessment
of Resistance of Medical Gloves
to Permeation by
Chemotherapy Drugs.

Meets and exceeds ASTM 6978 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs*

** Chemical permeation data of chemotherapy drugs and a Material Safety Data Sheet is available upon request.

Storage Recommendations:

Store in a cool dry place. Avoid excessive heat.

Open box should be stored away from direct sunlight, Fluorescent, lighting, moisture, x-rays and ozone conditions. Dispose in accordance with local regulations.

Quality Standards:

Manufactured in accordance with Quality System ISO 9001.

Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978.

This product is made from Synthetic Nitrile Polymer and contains no natural latex rubber.

* When used in accordance to the guidelines set forth in the USP 800 chapter.

Tested chemotherapy and Fentanyl drugs available upon request.



NitraCare 200

