



**MYCO**

**MEDICAL**

**TheSharpChoice<sup>®</sup>**

NMSDC Certified Diversity Supplier

**DATASHEET**

# SURGICAL N95 RESPIRATORS

NIOSH-approved Surgical N95 Respirators from MYCO Medical offer  $\geq 99.9\%$  bacterial and particle filtration efficiency\*, effective against particulate aerosols free of oil.

- 4 layers of protection
- FDA 510(k)-cleared; Class II medical device
- Meets ASTM F1862 Fluid Resistance Level 160 mmHg



- NIOSH Approved TC-84A-3348
- Cellulose-free; not made with natural rubber latex
- Flat folded disposable face mask with a pleated front and conformable nose clamp. This forms a cup-shape mask covering the wearers nose, mouth and chin.
- Constructed with non-woven fabric; orange on the outer side and white on the inner side.

- Elastic bands are ultrasonically welded and are placed at the crown and neck for comfort. No ear loops.
- Available in 2 sizes to fit a broad range of faces
- Country of origin: Japan; TAA Compliant



**Bringing Reliability, Innovation and Value to Healthcare<sup>®</sup>**



Front and side view of surgical N95 respirator being worn

## Surgical N95 Respirators

PRODUCT CODES	DESCRIPTION	UNIT OF MEASURE
MYCMYTN01-11	Medium N95 Respirator	50/bx, 6 bx/case
MYCMYTN01-12	Small N95 Respirator	50/bx, 6 bx/case

### List of FDA-cleared NIOSH-approved surgical N95 respirators:

[https://www.cdc.gov/niosh/npptl/topics/respirators/disp\\_part/resresource3surgicaln95.html](https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/resresource3surgicaln95.html)

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.