

DISCOVER THE CLEARSPEC DIFFERENCE

ClearSpec Disposable Speculum With Repposable LED Illuminator



THE ONLY ILLUMINATED SINGLE-USE VAGINAL SPECULUM TO FEATURE A REPOSABLE LED LIGHT SOURCE

ClearSpec's advanced design provides a clear view of the cervix and is ideal for use in many settings requiring vaginal exposure, access or biopsy. Additionally, the ClearSpec disposable speculum eliminates the need for reprocessing metal specula, and minimizes any risk of cross-contamination. The cordless LED illuminator provides crisp, natural white light designed for clinical ease of use.

CLEARSPEC BENEFITS



Reduces the risk of cross-contamination



Eliminates the expense of reprocessing metal specula



Offers savings when compared to other illuminated specula



Repposable LED illuminator can be used for up to 100 exams



Improves physician workflow




Quality patient experience

REPOSABLE ILLUMINATOR

Patented LED illuminator is designed to offer better balance, along with a crisp, natural white light for 2+ hours.

Cordless LED illuminator does not require recharging and automatically turns on when inserted into handle, and off when removed.

One FREE illuminator is included in each box of 24 specula.



A True Everyday Speculum



Ideal for all facilities and settings that require a vaginal speculum:

PELVIC +
CERVICAL EXAMS
ESSURE

HYSTEROSCOPY
CERVICAL
BIOPSIES

PAP SMEARS +
CERVICAL CULTURES
CRYOTHERAPY

ENDOMETRIAL
BIOPSY + ABLATION
COLPOSCOPY

ClearSpec Is the Clear Choice

BETTER PERFORMANCE AT A BETTER PRICE



The ClearSpec Disposable Speculum can offer significant savings compared to other disposable illuminated specula.



To add to the value, a free reusable LED illuminator is provided with each 24-count box and can be used for up to 100 exams.

AVAILABLE IN 3 CONVENIENT,
COLOR-CODED SIZES



ORDERING INFORMATION

CLEARSPEC DISPOSABLE VAGINAL SPECULUM

Graves Style, Non-sterile

102NS – Small (24 specula/box + 1 FREE illuminator, 96 specula/case)

103NS – Medium (24 specula/box + 1 FREE illuminator, 96 specula/case)

104NS – Large (18 specula/box, 72 specula/case)

CLEARSPEC REPOSABLE LED SPECULUM ILLUMINATOR

105 – Battery-operated LED light source (1/box, 12/case)