



Bemis Health Care develops innovative, thoughtfully-designed products to complement Hi-Flow™ suction canisters that are staff-friendly, efficient and safety-focused to save time and money.



Quick-Connect II

An all-in-one direct-to-regulator bracket solution designed for fast, safe and secure suction canister connections ... every time



Hi-Flow™ Suction Canister

High flow rate and suction dependability



Quick-Drain™ System

Designed to provide an efficient solution for control of safety and costs in the disposal of liquid infectious waste.



Bemis Fluid and Waste Management Solution

What is the Bemis Quick-Drain™ System?

An innovative, simple system for the safe disposal of infectious liquid waste that provides significant savings on solidifier/waste hauling costs.

It's as simple as 1-2-3 to save money on medical waste disposal



Call your local AMS Rep to learn how the Bemis Fluid and Waste Management Solution can save you money!

AtlanticMedicalSolutions.com
704-561-0004



Take filled suction canisters to Quick-Drain installed in existing utility room. System is low profile with a small wall space footprint.



NO electricity needed. System empties canisters using water pressure to create a simple vacuum. Requires only basic water/sanitary lines.



Engage rinse, then discard empty canisters – saves red bag waste hauling costs due to ~90% less weight / no solidifier – reduces solidifier and aerosolization exposure risks.