



The easiest way to test for lead

Test, educate and take action, all in one visit.



“ We switched to LeadCare® II and love it! We get results right away and don't have to contact patients for follow-up. ”


- Dr. Ben Spitalnick, Pediatrician

- **Convenient**
CLIA-waived, just two drops of blood.
- **Patient satisfaction**
Parents appreciate a single visit and a simple fingerstick.
- **Immediate action**
Results in 3 minutes, while the patient waits.
- **Saves time**
Eliminates the need for post-visit follow up.
- **Boosts compliance**
Ensures the test gets done.
- **Reimbursable**
CPT code 83655


What your peers love about leadcare® II

Ease of Use | Workflow Efficiency | Accuracy


TESTING IS AS EASY AS 1, 2, 3



1
Collect: 50 µL capillary whole blood (tubes included in the test kit).



2
Prepare: Dispense blood into a treatment reagent tube and mix.



3
Test: Place a drop of sample on the sensor; results in 3 minutes.

“ Increased screening significantly. Paid for itself in 9 months. ”


- Laura Rae Morette, CPNP

“ We love the test - quick, easy, accurate! ”

- Debra Mykkanen, FNP

ACCURATE

Compares to GFAAS Reference Method*



WE MAKE STARTUP SIMPLE

- Online video training & certification.
- Technical phone support for troubleshooting.
- Free patient educational materials.
- All you need is an analyzer & a test kit (includes materials & controls for 48 tests).

“ Easy set up, easy testing & fast results. ”

- Elide S. Silva, RN, BSN

LeadCare® II Analyzer Catalog Number: 70-6760

LeadCare® II Blood Lead Test Kit Catalog Number: 70-6762

*LeadCare® II Blood Lead Test Kit Package Insert