



ATLANTIC
MEDICAL
SOLUTIONS

**Modular, Multiplatform,
Affordable – Evolved**

AdView²

Vital Signs Monitors



*Some assembly required



www.adctoday.com/adview2



9005BP
BP Module

9005BPTO
BP/Temp

9005BPS
BP/SpO2

9005BPSTO
BP/SpO2/Temp

Our truly modular design allows you to create the diagnostic station that's perfectly suited to your needs. Start with the fully featured Blood Pressure base unit and add optional Temperature and SpO2 functionality, now or in the future. Order modules at the time of purchase and we'll ship your Adview 2 configured, or add them as you need them to upgrade in the field.

Ordering Information

Unit Models		Mounting Accessories and Platforms	
BP Only	9005BP	Desk Caddy*	9005D
BP/Temp	9005BPTO	Wall Mount	9005W
BP/SpO2	9005BPS	Mobile Stand	9006M
BP/Temp/SpO2	9005BPSTO	*Optional desktop caddy keeps cuffs and cords organized.	
Accessories	Barcode Scanner	WiFi Dongle	
Item No.	9005SCANL	9005DON	
	9005SCAN2D		

Customers ordering in Canada can call 1-800-ADC-2670 for model availability.



Included Adcuff+ Sizes:



Optional Adcuff+ Sizes:



Visit our website for all available accessories

Adview² Modular Diagnostic Station

The world's first — and still only — truly modular diagnostic station. Now in its second generation.

The Adview² features

- A flexible design with up to three modules: blood pressure, temperature, and pulse oximetry
- Clinical-grade Advantage™ BP technology from SunTech Medical®
- A proprietary manual BP mode, allowing measurement with a stethoscope using the auscultatory method
- New, customizable average mode, with setting options ranging from two to five readings and start and interval delays
- A temperature module featuring Covidien Filac™ technology, providing accurate measurement in about 10 seconds. Module comes with standard oral/axillary probe; optional rectal probe available.
- ChipOx Nellcor™ SpO2 module, providing accurate readings in about 10 seconds. Optional pediatric sensor available.
- Contemporary, compact design with an easy-to-read LCD display and an integrated handle for portability
- Lithium ion power for longer battery life and lighter weight
- Desktop, mobile, and wall mount options
- Optional accessories, including a Wi-Fi dongle and two scanner options (linear or 2D) for uploading patient data to your facility's EMR
- Network-enabled technology with IHE HL7 messaging
- Three-year warranty on the blood pressure base unit, SpO2 module, and temperature module; three-year warranty on cuffs; 90-day warranty on all other accessories



Adview 2



Shown here: Model 5610-39W
with Adview 2
(Model 9005BPT)

Diagnostix Adstation™ with Adview 2

3.5V Modular Diagnostix Wall System
with Adview 2

- Adview 2 modular diagnostic station with optional pulse oximetry and thermometry
- 3.5V Diagnostix Coax Plus Ophthalmoscope or Coax Ophthalmoscope
- 3.5V Diagnostix PMV Otoscope or Standard Otoscope with specula dispenser
- Add a 429 Non-Contact Thermometer
- Add a 3.5V Diagnostix Throat Illuminator or Dermoscope
- Choose from Halogen/Xenon or high-performance AdLED lighting technology

AtlanticMedicalSolutions.com

704-561-0004