



# SAFE life™



ATLANTIC  
MEDICAL  
SOLUTIONS

## T-Cup Multi-Drug Urine Test Cup

- CLIA Waived, 510(K) Cleared & Forensic Use Only
- Female User Friendly & Customizable
- Test up to 26 Drugs per Cup



### CLIA-waived Drugs & Cutoff Levels

	Amphetamine (AMP)	1000, 500 ng/mL		Tricyclic Antidepressants (TCA)	1000 ng/mL
	Barbiturates (BAR)	300 ng/mL		Methamphetamine (MET)/mAMP	1000, 500 ng/mL
	Benzodiazepines (BZO)	300 ng/mL		Methylenedioxymethamphetamine (MDMA)	500 ng/mL
	Buprenorphine (BUP)	10 ng/mL		Morphine (MOP/OPI300)	300 ng/mL
	Cocaine (COC)	300, 150 ng/mL		Opiate (OPI)	2000 ng/mL
	Cannabinoids (THC)	50 ng/mL		Oxycodone (OXY)	100 ng/mL
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL		Phencyclidine (PCP)	25 ng/mL
	Methadone (MTD)	300 ng/mL		Propoxyphene (PPX)	300 ng/mL

### Forensic Use Panels

	Monoacetylmorphine (6-MAM)	10 ng/mL		Synthetic Cannabinoid (K2)	50 ng/mL
	Amphetamine (AMP)	300, 500 ng/mL		Ketamine (KET)	300, 1000 ng/mL
	Buprenorphine (BUP)	5 ng/mL		Kratom (KRA)	100, 300 ng/mL
	Benzodiazepines (BZO)	100, 200 ng/mL		Lysergic Acid Diethylamide (LSD)	20 ng/mL
	Cocaine (COC)	100, 150 ng/mL		Methamphetamine (MET)/mAMP	300, 500 ng/mL
	Cotinine (COT)	200 ng/mL		Morphine (MOP or OPI)	100 ng/mL
	EDDP	100 ng/mL		Methadone (MTD)	200 ng/mL
	Ethyl Glucuronide (EtG)	500 ng/mL		Methaqualone (MQL)	300 ng/mL
	Fentanyl (FTY)	20 ng/mL		Pregabalin (PGB)	500 ng/mL
	Gabapentin (GAB)	2000 ng/mL		Cannabinoids (THC)	15, 25, 40 ng/mL
	Hydromorphone (HMO)	300 ng/mL		Tramadol (TRA)	100, 200, 1000 ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



# SAFE life™

## Product Information

※ Kit Size: 25 / Case Size: 4

Item Code	Description	Complexity
TDOA-154	5-Panel T-Cup; AMP, COC, OPI, PCP, THC	CLIA WAIVED
TDOA-254	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-254A3	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC with 3 AD	CLIA WAIVED
TDOA-164	6-Panel T-Cup; AMP, COC, MET, OPI, PCP, THC	CLIA WAIVED
TDOA-264	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-274	7-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
TDOA-3104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA WAIVED
TDOA-3104A3	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with 3 AD	CLIA WAIVED
TDOA-4104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-7104	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA WAIVED
TDOA-8104	10-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA WAIVED
TDOA-3124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-4124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, PCP, PPX, THC	CLIA WAIVED
TDOA-6124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-6124A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-6125A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125V	12-Panel T-Cup w/Vacu Lid; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-7125	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-7125A3	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-2135	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-2135A3	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-6135F	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA	Forensic Use Only
TDOA-7135A3EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, ETG, FTY, TRA with 3 AD	Forensic Use Only
TDOA-8135E	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA WAIVED
TDOA-1144A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-1145A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-8145EF	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only
TDOA-9145EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2	Forensic Use Only
TDOAC-L1223A6EFTKK6	22-Panel T-Cup Clear; AMP300, BAR300, BUP10, BZO200, COC100, EDDP300, HMO300, KET1000, MDMA500, MET300, MOP100, MTD200, OXY100, PCP25, PPX300, THC25, ETG500, FTY20, TRA200, K2 50, KRA300, 6-MAM10 +(CR, GL, NI, OX, pH, S.G.)	Forensic Use Only
FDOA-6125	12-Panel F-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED

Rel.: 2023/05/09