



Integrated Solution

OPTIMIZING SPACE IN THE ENT OFFICE

AMS

ORDER HERE • SALES REP

AMS

NEW ACCOUNT APPLICATION



OP-S9 Chair

OP-C16 Microscope

Otosmart Plus Treatment Cabinet

- Comfort and excellence for doctor and patient -
- High precision examinations and procedures -
- Compact and modern design -

OPTOMIC

producing quality

OTOSMART PLUS

Treatment Cabinet

OPTIMIZE SPACE AND IMPROVE WORK ERGONOMICS



New feature

The microscope can be mounted now on the OTOSMART PLUS cabinet

Professional - Quality imaging

- VESA Mount for Medical Monitor

Easy Disinfection

- Medical Grade Glass, obstacle-free surface for thorough cleaning

Integration of Endoscopes

- Integrated supports for rigid and flexible endoscopes, fiber optic cable and camera head

Versatility

Space for endoscopy equipment and medical devices

Mobility

- 4 omnidirectional wheels
- Mobility for cleaning and positioning for use

Integrated Suction and Pressure


- Automatic On/Off
- Suction: 70 L/min
- Pressure: 2 bar

Consultation Instruments

- 3 drawers for surgical instruments & supplies

Dicards

- Trash Bin
- Used instrument tray



 ORDER HERE • SALES REP

 NEW ACCOUNT APPLICATION

OP-S9

Procedures Chair

COMFORT AND EXCELLENCE FOR DOCTOR AND PATIENT



Tilt

- Visible precision
- 30° tilt memory (vestibular)
- Horizontal position
- Trendelenburg memory

Maximum load capacity

- 727 pounds
- Robust electric motor system

Electric brake

- Locks rotation with the press of a button

Exceptional Stability

- Sturdy construction with a firm and secure backrest
- Stable base with a compact design

Multiple handling options

- Side-mounted keyboard
- Pedal (optional)
- Handheld remote (optional)

Versatile Headrest

- Angle and height adjustable
- Quick-release system for easy vestibular testing
- Adjustable for children as well

Comfortable Cushions

- Range of 5 colors
- Pur Upholstery: plush and comfortable

Great rotation and smooth turning

- 380° rotation
- Agility for examining both ears
- Precise maneuvers without vibrations

Minimum seat height

- Minimum 19"
- Maximum 33.5"
- Easy access for patients with reduced mobility or disabilities

Seamless

- Facilitates cleaning
- Dirt does not penetrate the cushion

 **ORDER HERE • SALES REP**

 **NEW ACCOUNT APPLICATION**