



Bemis Health Care develops **innovative, thoughtfully-designed products** to complement **Hi-Flow™ suction canisters** that are **staff-friendly, efficient and safety-focused** to save time and money.



### Quick-Connect II



An **all-in-one direct-to-regulator** bracket solution designed for **fast, safe and secure suction canister connections ... every time.**

### Hi-Flow™ Suction Canister



High flow rate and suction dependability

### Quick-Drain™ System



Designed to provide an **efficient solution** for **control of safety and costs** in the disposal of liquid infectious waste.

**Bemis Fluid and Waste Management Solution**

## What is the Bemis Quick-Drain™ System?



Call your local AMS Rep to learn how the Bemis Fluid and Waste Management Solution can save you money!



An **innovative, simple system** for the **safe disposal of infectious liquid waste** that **provides significant savings** on solidifier/waste hauling costs.

It's as simple as 1-2-3 to save money on medical waste disposal



Take filled suction canisters to Quick-Drain installed in existing utility room. System is low profile with a small wall space footprint.



NO electricity needed. System empties canisters using water pressure to create a simple vacuum. Requires only basic water/sanitary lines.



Engage rinse, then discard empty canisters – saves red bag waste hauling costs due to ~90% less weight / no solidifier – reduces solidifier and aerosolization exposure risks.

Compatible with: Bemis Hi-Flow 1200cc, 2000cc 3000cc | Quick-Fit Liners 1000cc, 1500cc, 3000cc