



5 REASONS to Choose

semper med velvet™



velvet data sheet

3g Examination Nitrile | XS - XL (SMNP301-305) | 300/bx (250XL)

2

**LIGHTWEIGHT,
BUT STRONG
3 GRAM**



As a true manufacturer, we control the formula of the glove. Nitrile is a superior glove material boasting excellent chemical resistance, dexterity, and tactile sensitivity.



3

**BETTER VALUE
300 GLOVES /BOX BY WEIGHT
(250 XL)**



300 gloves per box means greater storage efficiencies and less packaging waste.

1

**MADE IN
MALAYSIA**

We own two Malaysian factories utilizing state-of-the-art technology to achieve the highest quality in our medical and general purpose grade disposable gloves.

Central Medicare
Teluk Intan, Malaysia



4

**SOFT AS VELVET
ENHANCED
COMFORT**

The ultra-soft, ultra-light nitrile innovative design results in a glove that conforms to the hand and offers a comfortable, natural fit with less hand fatigue and unrestricted hand motion.

5

**TESTED AGAINST
56 CHEMO DRUGS*,
FENTANYL CITRATE,
& GASTRIC ACID**



Tested for use with Chemotherapy drugs and Fentanyl Citrate as per ASTM D6978. Meets USP 800 requirements for PPE compliance for chemotherapy.

* Testing results available by request.
Not recommended for handling Carmustine or Thiotepa

FREE Samples!
Contact your local representative



rev02072025

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is available at www.gotgloves.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. HARPS is not liable for incorrect use of the gloves nor will HARPS be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used. In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. HARPS bears no liability for allergic reactions. ©Copyright 2024 HARPS

