

OUR COMMITMENT TO QUALITY, RESILIENCY AND EXCELLENCE

For over 50 years, Bemis Health Care has been known for providing quality medical device products, excellent customer service and supply resiliency. We are proudly part of Bemis Manufacturing Company, a fourth generation family-owned US business founded in 1901. Bemis Health Care recognizes that while customers have many choices in the products they use, we strongly believe what sets us apart from the rest is our ability to provide:

QUALITY PRODUCTS — Bemis Health Care is committed to and takes great pride in manufacturing quality products. We are dedicated to providing safe, dependable and thoughtfully-designed products that are clinician friendly, efficient and safety focused.

OUTSTANDING CUSTOMER SERVICE — Customer support at Bemis Health Care is a full team effort. We work together to address your needs in a timely, professional and responsive manner. Our customer relationships are true partnerships and we see challenges as opportunities for Bemis Health Care to make the difference in your experience.

SUPPLY RESILIENCY — In today's healthcare landscape, working with reliable and adaptable suppliers is essential for maintaining a healthy supply chain. Bemis Health Care can quickly adjust to changing market conditions, fluctuating demand and customer expectations because of our locally-sourced materials and American manufacturing capabilities.



[BemisHealthCare.com](https://www.BemisHealthCare.com)



BEMIS HEALTH CARE PROVIDES HIGH-QUALITY MEDICAL PRODUCTS FOR THE HEALTHCARE INDUSTRY

Bemis offers high quality, single-use containers available in a variety of sizes and clinical use types for biohazardous sharps, trace chemotherapy, non-hazardous and hazardous pharmaceutical waste disposal. Accessory solutions available to fit a variety of mounting and mobility needs.



BIOHAZARD SHARPS CONTAINERS

Bemis Sharps Containers are designed to meet OSHA requirements for safe collection of biohazardous sharps. Containers are available for use in both public access and restricted clinical areas. Biohazard logo clearly distinguishes the container use from other disposal collection types. Available in multiple sizes, colors and lid openings to meet clinical use needs.



CHEMOTHERAPY SHARPS CONTAINERS

Bemis Chemotherapy Containers are designed for safe collection and containment of trace chemotherapy waste. Trace waste is defined by CDC/NIOSH as "containing less than 3% by weight of original quantity of hazardous drug." Each container includes Biohazard and Cytotoxic logos indicating clear disposal type. All three sizes are available in both yellow and white colors and include gasketed lids with absorbent pads.



NON-HAZARDOUS PHARMACEUTICAL CONTAINERS




















Bemis Non-Hazardous Pharmaceutical Containers are designed for the disposal of expired or unused non-hazardous pharmaceuticals. Each container includes gasketed lids with absorbent pads. Available in three blue container sizes with labeling to easily distinguish from other disposal collection types.



HAZARDOUS PHARMACEUTICAL CONTAINERS

Bemis Hazardous Pharmaceutical Waste Containers are designed for the collection of P and U-listed chemicals according to the EPA Resource Conservation and Recovery Act (RCRA). Each container includes gasketed lids with absorbent pads. Available in three black container sizes with labeling to easily distinguish from other disposal collection types.

BEMIS HEALTH CARE - WASTE DISPOSAL SOLUTIONS

TYPE	CAPACITY	CATALOG#	COLOR	LID TYPE	SAMPLE WASTE
BIOHAZARD SHARPS*					
Phlebotomy Sharps	1-Quart	100 030	 Red	Dual openings/ vertical drop	1-10 mL syringes, blood collection tubes, pipettes
Patient/Exam Room Sharps	5-Quart	175 030	 Red	Rotating cylinder/ horizontal drop	1-60 mL syringes with needle
		175 020	 Beige		
Multi-Purpose Sharps	2-Gallon	102 030	 Red	Angled chute/ horizontal drop	1-60 mL syringes with needle, blood collection tubes, slides, pipettes
		102 020	 Beige		
Patient/Exam Room Sharps	3-Gallon	333 030	 Red	Rotating cylinder/ horizontal drop	1-60 mL syringes with needle
		333 020	 Beige		
Multi-Purpose Sharps	3 Gallon	303 030	 Red	Angled chute/ horizontal drop	1-60 mL syringes with needle, blood collection tubes, slides, pipettes
		303 020	 Beige		
Multi-Purpose Large Volume Sharps	8-Gallon	108 030	 Red	Hinged for large items/ separate opening for smaller items	Various sharps and larger OPIM items
	11-Gallon	111 030	 Red		
TRACE CHEMOTHERAPY SHARPS					
Chemotherapy Medical Waste	2-Gallon	202 004	 White	Screw-top with 4-inch opening	Various trace chemo sharps, vials, syringes
		202 040	 Yellow		
Chemotherapy Large Volume Medical Waste	8-Gallon	208 004	 White	Gasketed hinge for large items/separate opening for smaller items	Various trace chemo sharps, tubing, vials, syringes, pads, gowns
		208 040	 Yellow		
	11-Gallon	211 004	 White		
		211 040	 Yellow		
NON-HAZARDOUS PHARMACEUTICAL					
Pharmaceutical Medical Waste	2-Gallon	4002 050	 Blue	Screw-top with 4-inch opening	NO SHARPS - Non-hazardous pharmaceutical syringes/vials with residual medication, small items
Pharmaceutical Large Volume Medical Waste	8-Gallon	4008 050	 Blue	Hinged for large items/ separate opening for smaller items	NO SHARPS - Non-hazardous pharmaceutical syringes/vials with residual medication, tablets, creams, eye drops
	11-Gallon	4011 050	 Blue		
HAZARDOUS PHARMACEUTICAL					
Hazardous Pharmaceutical Medical Waste	2-Gallon	5002 070	 Black	Screw-top with 4-inch opening	NO SHARPS - Small hazardous pharmaceutical P and U-listed items
Hazardous Pharmaceutical Medical Waste	8-Gallon	5008 070	 Black	Hinged for large items/ separate opening for smaller items	NO SHARPS - Hazardous pharmaceutical P and U-listed items
	11-Gallon	5011 070	 Black		

* Yellow Biohazard Sharps Containers for Canadian market not listed in table - refer to pages 6-8 for additional information

GENERAL MEDICAL WASTE DISPOSAL + SEGREGATION GUIDE

CHECK STATE OR LOCAL GUIDELINES FOR SPECIFIC MEDICAL WASTE DISPOSAL REGULATIONS

LANDFILL



BAG WASTE - WHITE OR CLEAR

- IV bags and tubing
- Empty med vials
- Trash/wrappers
- Dressings, Chux
- Empty drainage bags
- Disposable items
- Sanitary napkins
- Batteries

AUTOClave / INCINERATION



RED BAG WASTE - BIOHAZARD

- Infectious materials (OPIM)
- Blood tubing, bags, Hemovac®, Pleur-evac®
- Blood soaked or dripping dressings
- Contaminated waste from isolation patients
- Suction canisters and liners with bodily fluids



SHARPS - BIOHAZARD

- Needles
- Glass slides, pipettes
- Scalpels
- Bloody broken glass
- Used sutures
- Empty syringes
- Trocars, introducers, catheters, electrocautery tips

LOW-TEMP INCINERATION



TRACE CHEMO - BIOHAZARD

- Trace chemo waste of no more than 3% weight of original components
- All supplies used to make and administer chemotherapy medication including sharps
- Tubing, empty bags/bottles/vials, syringes, needles, gloves, pads, masks, gowns, wipes



PHARMA - NON-HAZARDOUS

- NO SHARPS
- Residual (pourable) meds
- Discontinued IV bags/tubing
- Partially-used prescription or OTC medication
- Vials, tablets, capsules, powders, contained liquids, creams, lotions, eye drops, residual or wasted narcotics, unused narcotic patches

HIGH-TEMP INCINERATION



PHARMA - HAZARDOUS

- NO SHARPS
- Hazardous chemicals as defined by Resource Conservation and Recovery Act (RCRA) for P and U-listed chemicals
- Check the EPA website for current listing of P and U-listed hazardous waste

PATIENT + EXAM ROOM SHARPS CONTAINERS








DESIGNED FOR THE DISPOSAL OF BIOHAZARD SHARPS IN PATIENT AND EXAM ROOMS



FEATURES + BENEFITS

- Patient and Exam Room Sharps Containers are designed with a rotating lid that restricts access to contents and are used in patient or exam rooms, and other high traffic or public areas. Available in 5-Quart and 3-Gallon sizes with convenient mounting options, including enclosed locking cabinets, lockable wall brackets or tabletop use.
- Translucent containers and lids help identify fill level
 - Horizontal drop allows maximum capacity of sharps
 - Lid design accommodates syringes up to 60cc
 - Containers lock for final disposal
 - Nestable for space saving storage
 - English, Spanish and French labeling
 - Requires autoclave to be rendered non-infectious, then incinerate or grind/shred disposal

PATIENT + EXAM ROOM SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 5-QUART SHARPS CONTAINER <ul style="list-style-type: none">• Rotating cylinder top allows for hands-free disposal• Built-in handles• Uses Cabinet #1051 or Bracket/Key #445 for mounting HxWxL: 10 x 5 ¼ x 11 in	<div><div>■ 175 020</div>Translucent Beige</div> <div><div>■ 175 030</div>Translucent Red</div> <div><div>■ 175 040*</div>Translucent Yellow</div>	32/CASE	<div> 1051</div> <div> 445 020</div>
 3-GALLON SHARPS CONTAINER <ul style="list-style-type: none">• Rotating cylinder top allows for hands-free disposal• Built-in handles• Uses Cabinet #1031 or Bracket 423 or 425 for mounting HxWxL: 15 x 6 7⁄8 x 13 7⁄8 in	<div><div>■ 333 020</div>Translucent Beige</div> <div><div>■ 333 030</div>Translucent Red</div> <div><div>■ 333 040*</div>Translucent Yellow</div>	12/CASE	<div> 1031</div> <div> 423 000</div> <div> 425 000</div>

* Available for Canadian market

PHLEBOTOMY + MULTI-PURPOSE SHARPS CONTAINERS










DESIGNED FOR THE DISPOSAL OF BIOHAZARD SHARPS IN RESTRICTED AREAS



FEATURES + BENEFITS








- Multi-Purpose Sharps Containers are designed for use in restricted areas such as nursing stations, laboratories, operating rooms or pharmacies and available in multiple sizes to accommodate clinical needs.
- Built-in handles for easy transport (excluding phlebotomy)
 - Horizontal drop on most units maximizes container volume
 - Lids feature temporary closure or locked for disposal
 - Translucent containers and lids help identify fill level on 1-quart, 2 and 3-gallon sizes
 - 2 and 3-Gallon can be wall-mounted or used freestanding on tabletop or cart
 - 8 and 11-Gallon have large, hinged lid openings
 - Nestable for space saving storage
 - English, Spanish and French labeling
 - Requires autoclave to be rendered non-infectious, then incinerate or grind/shred disposal

PHLEBOTOMY + MULTI-PURPOSE SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 1-QUART PHLEBOTOMY CONTAINER <ul style="list-style-type: none">• Dual openings, vertical drop• Uses Holder #400 to secure on tray or tabletop HxWxL: 6 7/16 x 4 1/4 x 4 1/4 in	 100 030 Translucent Red  100 040* Translucent Yellow	100/CASE	 400 000
 2-GALLON SHARPS CONTAINER <ul style="list-style-type: none">• Accommodates syringes up to 60cc• Built-in handles• Uses Bracket/Key #445 for wall mounting HxWxL: 8 5/8 x 7 3/4 x 11 5/8 in	 102 020 Translucent Beige  102 030 Translucent Red  102 040* Translucent Yellow	30/CASE	 445 020

* Available for Canadian market

PHLEBOTOMY + MULTI-PURPOSE SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 <p>3-GALLON SHARPS CONTAINER</p> <ul style="list-style-type: none"> Accommodates syringes up to 60cc Built-in handles Uses Bracket #423 or 425 for wall mounting <p>HxWxL: 13 ½ x 6 ⅞ x 13 ⅞ in</p>	<p>■ 303 020 Translucent Beige</p> <p>■ 303 030 Translucent Red</p> <p>■ 303 040* Translucent Yellow</p>	<p>12/CASE</p>	 423 000  425 000
 <p>8-GALLON SHARPS CONTAINER</p> <ul style="list-style-type: none"> Hinged lid accommodates large items Separate round opening accepts small items Built-in handles Uses Cart #430 for mobility <p>HxWxL: 15 ⅞ x 11 ⅜ x 16 ½ in</p>	<p>■ 108 030 Red</p> <p>■ 108 040* Yellow</p>	<p>10/CASE</p>	 430 030
 <p>11-GALLON SHARPS CONTAINER</p> <ul style="list-style-type: none"> Hinged lid accommodates large items Separate round opening accepts small items Built-in handles Uses Cart #430 for mobility <p>HxWxL: 22 ½ x 11 ⅜ x 16 ½ in</p>	<p>■ 111 030 Red</p> <p>■ 111 040* Yellow</p>	<p>6/CASE</p>	 430 030

* Available for Canadian market

CHEMOTHERAPY SHARPS CONTAINERS







DESIGNED FOR THE DISPOSAL OF TRACE CHEMOTHERAPY WASTE



FEATURES + BENEFITS

- Designed for safe collection and containment of trace chemotherapy sharps and waste including needles, syringes, IVs, empty or partially empty vials, as well as gloves and gowns. Available in white or yellow and three sizes, Chemotherapy Containers are clearly labeled with cytotoxic and biohazard symbols to distinguish from other disposal types.
- Gasketed lids and absorbent pads provide leak resistance
 - Lids feature temporary closure or locked for disposal
 - 2-Gallon features a 4-inch screw top opening
 - 8 and 11-Gallon have large, hinged lid openings
 - Nestable for space saving storage
 - English, Spanish and French labeling
 - Requires low-temp incineration disposal

CHEMOTHERAPY SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 <p>2-GALLON CHEMOTHERAPY CONTAINER</p> <ul style="list-style-type: none">• Screw-on cap• Posts on lid for holding cap in place when in use• For use on tabletop or wall-mounted• Uses Bracket/Key #445 for wall mounting <p>HxWxL: 9 x 7 ¾ x 11 ⅝ in</p>	<div><input type="checkbox"/> 202 004 White w/gasket lid, absorbent pad</div> <div><input checked="" type="checkbox"/> 202 040 Yellow w/gasket lid, absorbent pad</div>	30/CASE	 <p>445 020</p>
 <p>8-GALLON CHEMOTHERAPY CONTAINER</p> <ul style="list-style-type: none">• Hinged lid accommodates large items• Separate round opening accepts small items• Built-in handles• Uses cart #430 for mobility <p>HxWxL: 15 7/8 x 11 13/16 x 16 1/2 in</p>	<div><input type="checkbox"/> 208 004 White w/gasket lid, absorbent pad</div> <div><input checked="" type="checkbox"/> 208 040 Yellow w/gasket lid, absorbent pad</div>	10/CASE	 <p>430 030</p>
 <p>11-GALLON CHEMOTHERAPY CONTAINER</p> <ul style="list-style-type: none">• Hinged lid accommodates large items• Separate round opening accepts small items• Built-in handles• Uses cart #430 for mobility <p>HxWxL: 22 1/2 x 11 13/16 x 16 1/2 in</p>	<div><input type="checkbox"/> 211 004 White w/gasket lid, absorbent pad</div> <div><input checked="" type="checkbox"/> 211 040 Yellow w/gasket lid, absorbent pad</div>	6/CASE	 <p>430 030</p>

NON-HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS

DESIGNED FOR THE DISPOSAL OF NON-HAZARDOUS PHARMACEUTICALS












FEATURES + BENEFITS

Pharmaceutical Waste Containers are designed for the disposal of expired or unused nonhazardous pharmaceuticals such as antibiotics, IVs, pills, etc. and are often used in pharmacies, nursing stations and med carts. The blue Pharmaceutical Waste containers include special labeling to easily distinguish from other disposal containers and are not for Sharps use.

- Gasketed lids and absorbent pads provide leak resistance
- Lids feature temporary closure or locked for disposal
- 2-Gallon features a 4-inch screw top opening
- 8 and 11-Gallon have large, hinged lid openings
- Nestable for space saving storage
- English, Spanish and French labeling
- Container NOT for Sharps use
- Requires low- temp incineration disposal

NON-HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 <p>2-GALLON PHARMACEUTICAL CONTAINER</p> <ul style="list-style-type: none"> • Screw-on cap • Posts on lid for holding cap in place when in use • For use on tabletop or wall-mounted • Uses Bracket/Key #445 for wall mounting <p>HxWxL: 9 x 7 ³/₄ x 11 ⁵/₈ in</p>	<p> 4002 050 Blue w/gasket lid, absorbent pad</p>	30/CASE	 <p>445 020</p>
 <p>8-GALLON PHARMACEUTICAL CONTAINER</p> <ul style="list-style-type: none"> • Hinged lid accommodates large items • Separate round opening accepts small items • Built-in handles • Uses Cart #430 for mobility <p>HxWxL: 15 ⁷/₈ x 11 ¹³/₁₆ x 16 ¹/₂ in</p>	<p> 4008 050 Blue w/gasket lid, absorbent pad</p>	10/CASE	 <p>430 030</p>
 <p>11-GALLON PHARMACEUTICAL CONTAINER</p> <ul style="list-style-type: none"> • Hinged lid accommodates large items • Separate round opening accepts small items • Built-in handles • Uses Cart #430 for mobility <p>HxWxL: 22 ¹/₂ x 11 ¹³/₁₆ x 16 ¹/₂ in</p>	<p> 4011 050 Blue w/gasket lid, absorbent pad</p>	6/CASE	 <p>430 030</p>

HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS

DESIGNED FOR THE DISPOSAL OF HAZARDOUS PHARMACEUTICALS







FEATURES + BENEFITS

Hazardous Pharmaceutical Waste Containers are designed to meet EPA requirements of disposing P and U-listed chemicals. Black containers are available in three sizes, include special labeling to easily distinguish from other disposal containers and are not for Sharps use.







- Gasketed lids and absorbent pads provide leak resistance
- Lids feature temporary closure or locked for disposal
- 2-Gallon features a 4-inch screw top opening
- 8 and 11-Gallon have large, hinged lid openings
- Nestable for space saving storage
- English, Spanish and French labeling
- Container NOT for Sharps use
- Requires high-temp incineration disposal

HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING	COMPATIBLE
 <p>2-GALLON HAZARDOUS PHARMACEUTICAL</p> <ul style="list-style-type: none"> • Screw-on cap • Posts on lid for holding cap in place when in use • For use on tabletop or wall-mounted • Uses Bracket/Key #445 for wall mounting <p>HxWxL: 9 x 7 ³/₄ x 11 ⁵/₈ in</p>	<p>■ 5002 070</p> <p>Black w/gasket lid, absorbent pad</p>	30/CASE	 <p>445 020</p>
 <p>8-GALLON HAZARDOUS PHARMACEUTICAL</p> <ul style="list-style-type: none"> • Hinged lid accommodates large items • Separate round opening accepts small items • Built-in handles • Uses Cart #430 for mobility <p>HxWxL: 15 ⁷/₈ x 11 ¹³/₁₆ x 16 ¹/₂ in</p>	<p>■ 5008 070</p> <p>Black w/gasket lid, absorbent pad</p>	10/CASE	 <p>430 030</p>
 <p>11-GALLON HAZARDOUS PHARMACEUTICAL</p> <ul style="list-style-type: none"> • Hinged lid accommodates large items • Separate round opening accepts small items • Built-in handles • Uses Cart #430 for mobility <p>HxWxL: 22 ¹/₂ x 11 ¹³/₁₆ x 16 ¹/₂ in</p>	<p>■ 5011 070</p> <p>Black w/gasket lid, absorbent pad</p>	6/CASE	 <p>430 030</p>

SHARPS CABINETS, MOUNTS + ACCESSORIES

CABINETS + ACCESSORIES

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING
 5-QUART WALL CABINET <ul style="list-style-type: none"> Two universal keys provided Easy mount "Swiss Cheese" style backing For quick replacement of Non-Bemis cabinets For use with sharps container #175 HxWxL: 11 ½ x 6 x 13 in 	1051	1/CASE
 3-GALLON WALL CABINET <ul style="list-style-type: none"> Two universal keys provided Easy mount "Swiss Cheese" style backing For quick replacement of Non-Bemis cabinets For use with sharps container #333 HxWxL: 16 ¾ x 7 x 15 in 	1031	1/CASE
 GLOVE BOX HOLDER <ul style="list-style-type: none"> Conveniently attaches to bottom of cabinet For use with wall cabinets #1051 and 1031 One glove box per cabinet HxWxL: 5 ½ x 4 ⅞ x 11 ⅝ in 	1041	1/CASE

CARTS + HOLDERS



ROLLING CART FOR 8 + 11-GALLON CONTAINERS

- Durable construction
- Accommodates all 8 and 11-gallon Containers

HxWxL: 16 x 13 x 17 in

430 030

1/CASE



PHLEBOTOMY HOLDER







- Translucent white allows visibility of container label
- Adhesive hook-and-loop secures holder in place
- For use with Phlebotomy Container #100

HxWxL: 6 ⅝ x 4 ¼ x 4 ¼ in

400 000

20/CASE

BRACKETS + ACCESSORIES

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING
 MOUNTING BRACKET/KEY FOR 5-QUART + 2-GALLON CONTAINERS <ul style="list-style-type: none"> • Mounting hardware included with 10 keys • Compatible with all Bemis 5-Quart and 2-Gallon Containers HxWxL: 6 x 2 ¼ x 10 ¼ in	445 020	10/CASE
 MOUNTING BRACKET FOR 5-QUART + 2-GALLON CONTAINERS <ul style="list-style-type: none"> • Mounting hardware included • Compatible with all Bemis 5-Quart and 2-Gallon Containers • Requires Key #440 to remove container from bracket HxWxL: 6 x 2 ¼ x 10 ¼ in	435 020	10/CASE
 PLASTIC KEY FOR #435 MOUNTING BRACKET <ul style="list-style-type: none"> • Unlocks 5-Quart and 2-Gallon Containers from Bracket #435 and 445 HxWxL: 3 ⅞ x ⅝ x 1 in	440 020	10/CASE
 NON-LOCKING WIRE BRACKET FOR 3-GALLON <ul style="list-style-type: none"> • Compatible with all Bemis 3-Gallon Containers • Mounting hardware included HxWxL: 7 ½ x 6 x 12 in	423 000	5/CASE
 LOCKING WIRE BRACKET FOR 3-GALLON <ul style="list-style-type: none"> • Compatible with all Bemis 3-Gallon Containers • With swing arm that fits over container and secures with padlock • Mounting hardware included • Requires Padlock/Key #426 HxWxL: 9 ½ x 6 x 12 in	425 000	5/CASE
 PADLOCK/KEY <ul style="list-style-type: none"> • For use with Locking Wire Bracket #425 HxWxL: 1 ¾ x ½ x 1 ⅞ in	426 000	1/CASE

QUICK REFERENCE TABLE - CONTAINER MOUNTING OPTIONS

			TABLETOP	MOUNT					CART
DESCRIPTION			PHLEBOTOMY	5-QUART	3-GALLON	WIRE	WIRE	WALL	ROLLING
MOUNT TYPE			CONTAINER HOLDER	SHARPS CABINET	SHARPS CABINET	NON-LOCKING BRACKET	LOCKING BRACKET*	KEY/LOCKING BRACKET**	SHARPS CART
PRODUCT#			400 000	1051	1031	423 000	425 000	445 020	430 030
BIOHAZARD SHARPS CONTAINERS									
100 040	1-Quart	Yellow	●						
100 030	1-Quart	Red	●						
175 040	5-Quart	Yellow		●				●	
175 030	5-Quart	Red		●				●	
175 020	5-Quart	Beige		●				●	
102 040	2-Gallon	Yellow						●	
102 030	2-Gallon	Red						●	
102 020	2-Gallon	Beige						●	
303 040	3-Gallon	Yellow				●	●		
303 030	3-Gallon	Red				●	●		
303 020	3-Gallon	Beige				●	●		
333 040	3-Gallon	Yellow			●	●	●		
333 030	3-Gallon	Red			●	●	●		
333 020	3-Gallon	Beige			●	●	●		
108 030	8-Gallon	Red							●
108 040	8-Gallon	Yellow							●
111 030	11-Gallon	Red							●
111 040	11-Gallon	Yellow							●
TRACE CHEMOTHERAPY SHARPS CONTAINERS									
202 004	2-Gallon	White						●	
202 040	2-Gallon	Yellow						●	
208 004	8-Gallon	White							●
208 040	8-Gallon	Yellow							●
211 004	11-Gallon	White							●
211 040	11-Gallon	Yellow							●
NON-HAZARDOUS PHARMACEUTICAL CONTAINERS									
4002 050	2-Gallon	Blue						●	
4008 050	8-Gallon	Blue							●
4011 050	11-Gallon	Blue							●
HAZARDOUS PHARMACEUTICAL CONTAINERS									
5002 070	2-Gallon	Black						●	
5008 070	8-Gallon	Black							●
5011 070	11-Gallon	Black							●

* Purchase Padlock and Key separately #426 000

** Wall Brackets and Keys can be purchased separately #435 020 and 440 020

NOTES

This image shows a full page of white paper with horizontal dotted lines. The lines are evenly spaced and run across the width of the page, providing a guide for handwriting practice. There are no margins, text, or other markings on the page.

