


S semper
med



ATLANTIC
MEDICAL
SOLUTIONS

roseTM

coated with aloe



**Coated with aloe to help
moisturize hands**

*Manufactured (tariff-free) in
Malaysia using state-of-the-art
technology to achieve quality that
exceeds FDA standards.*



Hospital



Beauty,
Body Art, &
Hairdressing



Dental



Primary
Care



Veterinary



HARPS
GLOBAL



www.gotgloves.com | 800.366.9545

roseTM coated with aloe

Sempermed® roseTM exam nitrile gloves are designed to elevate comfort and protection in any setting, from exam rooms to tattoo shops to day spa. These gloves feature a soothing aloe coating that helps prevent skin irritation, keeping hands moisturized and comfortable during extended use. The unique formulation of these gloves enhances protection while maintaining exceptional tactile sensitivity, allowing for precise handling without compromising comfort.

- **Coated with aloe:** Helps moisturize dry hands and helps prevent skin irritation while you work. A great choice for sensitive skin.
- **Fun, unique color:** Unique rose pigment color brightens any environment from the ER to the day spa.

- **Protects users and patients from the transmission of infections:** Designed to provide protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

Technical data

TYPE	Nitrile single use disposable protective glove, non-sterile	
GRADE	Examination	
COLOR	Fuchsia	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Finger textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	3 years in original packaging if stored as per DIN 7716, ISO 2230	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger .10mm 3.9mils	Palm .06mm 2.4mils based on average production value measured in single layer
TENSILE STRENGTH (after aging)	14MPa (min)	
ELONGATION (after aging)	500% (min)	
BARRIER PERFORMANCE	AQL ≤ 1.5	
TESTING	ASTM D412, ASTM D573, ASTM D3767, ASTM D5151 and ASTM D6124	

Purpose

Single use medical device class 1, as per the FDA Single-use medical device

Suitable for food handling, as per relevant sections of 21 CFR Parts 170-199

Packaging

DIMENSIONS	Dispenser 245L x 125W x 74H mm Carton 385L x 255W x 255H mm
------------	--

Reorder numbers / Glove counts*

SMNR201	X-Small	200 gloves/bx	10 bxs/cs
SMNR202	Small	200 gloves/bx	10 bxs/cs
SMNR203	Medium	200 gloves/bx	10 bxs/cs
SMNR204	Large	200 gloves/bx	10 bxs/cs
SMNR205	X-Large	180 gloves/bx	10 bxs/cs

SAMPLES: Available by request - ask your sales rep!

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes

*glove counts "by weight"



AtlanticMedicalSolutions.com
704-561-0004

rev 04/2025 US

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is available at www.gotgloves.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. HARPS is not liable for incorrect use of the gloves nor will HARPS be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used. In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. HARPS bears no liability for allergic reactions. ©Copyright 2024 HARPS